

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/13/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 06/12/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 06/13/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131205

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2851-01	KQA067F-1-1	1	1.00	1.00	2.00	2.00	100.0	oilc
P2851-02	KQA067F-1-2	2	1.00	1.00	2.00	2.00	100.0	oilc
P2851-03	BC278423-1-1	3	1.00	1.00	2.00	2.00	100.0	oilc
P2851-04	BC278423-1-2	4	1.00	1.00	2.00	2.00	100.0	oilc
P2851-05	BC247727-1-1	5	1.00	1.00	2.00	2.00	100.0	oilc
P2851-06	BC247727-1-2	6	1.00	1.00	2.00	2.00	100.0	oilc
P2852-01	31553	7	1.00	1.00	2.00	2.00	100.0	wipe sample
P2852-02	31554	8	1.00	1.00	2.00	2.00	100.0	wipe sample
P2852-03	31555	9	1.00	1.00	2.00	2.00	100.0	wipe sample
P2852-04	31556	10	1.00	1.00	2.00	2.00	100.0	wipe sample
P2852-05	31557	11	1.00	1.00	2.00	2.00	100.0	wipe sample
P2853-01	1	12	1.18	8.65	9.83	9.15	92.1	
P2853-02	2	13	1.17	8.54	9.71	9.24	94.5	
P2853-03	3	14	1.18	8.45	9.63	9.08	93.5	
P2853-04	4	15	1.16	8.75	9.91	9.43	94.5	
P2853-05	5	16	1.13	8.81	9.94	9.9	99.5	
P2853-06	6	17	1.16	8.58	9.74	7.99	79.6	
P2853-07	7	18	1.16	8.47	9.63	8.22	83.4	
P2853-08	8	19	1.14	8.58	9.72	9.31	95.2	
P2853-09	9	20	1.16	8.83	9.99	9.47	94.1	
P2853-10	10	21	1.13	8.75	9.88	8.8	87.7	
P2853-11	11	22	1.15	8.62	9.77	8.99	91.0	
P2853-12	12	23	1.14	8.45	9.59	8.92	92.1	
P2853-13	13	24	1.16	8.83	9.99	9.66	96.3	
P2855-01	BC278406-1-1	25	1.00	1.00	2.00	2.00	100.0	oilc
P2855-02	BC278406-1-2	26	1.00	1.00	2.00	2.00	100.0	oilc
P2855-03	BC274767-1-1	27	1.00	1.00	2.00	2.00	100.0	oilc
P2855-04	BC274767-1-2	28	1.00	1.00	2.00	2.00	100.0	oilc

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/13/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 06/12/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 06/13/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131205

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2855-05	72-11989	29	1.17	8.50	9.67	8.65	88.0	
P2855-06	72-11989-E2	30	1.15	8.51	9.66	7.99	80.4	
P2855-07	ARIEL-CABLE-1-1	31	1.00	1.00	2.00	2.00	100.0	oilc
P2855-08	ARIEL-CABLE-1-2	32	1.00	1.00	2.00	2.00	100.0	oilc
P2855-09	ARIEL-CABLE-1-3	33	1.00	1.00	2.00	2.00	100.0	oilc
P2855-10	ARIEL-CABLE-1-4	34	1.00	1.00	2.00	2.00	100.0	oilc
P2855-11	ARIEL-CABLE-1-5	35	1.00	1.00	2.00	2.00	100.0	oilc
P2855-12	ARIEL-CABLE-1-6	36	1.00	1.00	2.00	2.00	100.0	oilc
P2857-01	72-11985	43	1.00	1.00	2.00	2.00	100.0	concrete sample, 100% solids
P2857-02	ETGI-355	44	1.00	1.00	2.00	2.00	100.0	concrete sample, 100% solids
P2858-01	NB-07-06122024	45	1.13	8.51	9.64	8.96	92.0	
P2858-02	NB-07-06122024-E2	46	1.12	8.82	9.94	9.18	91.4	
P2859-01	OK-02-06122024	47	1.13	8.41	9.54	8.9	92.4	
P2859-02	OK-02-06122024-E2	48	1.12	8.62	9.74	8.88	90.0	
P2864-01	WC-H	37	1.17	8.41	9.58	8.37	85.6	
P2864-02	WC-H-EPH	38	1.18	8.71	9.89	8.5	84.0	
P2864-03	WC-H-VOC	39	1.19	8.59	9.78	8.9	89.8	
P2864-05	WC-I	40	1.19	8.61	9.8	8.18	81.2	
P2864-06	WC-I-EPH	41	1.2	8.69	9.89	8.79	87.3	
P2864-07	WC-I-VOC	42	1.2	8.58	9.78	9.11	92.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$