

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/25/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 06/24/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/25/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB131387

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3040-01	DA053P-SD02-062124-00	1	1.18	8.49	9.67	7.77	77.6	
P3040-02	DA053P-SD01-062124-00	2	1.13	8.80	9.93	8.5	83.7	
P3040-03	SS050P-SD03-062124-00	3	1.15	8.80	9.95	8.39	82.3	
P3040-04	SS050P-SD02-062124-00	4	1.15	8.82	9.97	7.79	75.3	
P3040-05	SS030P-SD01-062124-00	5	1.18	8.40	9.58	8.01	81.3	
P3040-06	ST009P-SD03-062124-00	6	1.16	8.75	9.91	5.72	52.1	
P3041-01	FT013P-SD01-062024-00	7	1.16	8.41	9.57	7.56	76.1	
P3041-02	SS030P-SD03-062024-00	8	1.16	8.57	9.73	5.21	47.3	
P3041-03	FT013P-SD02-062024-00	9	1.17	8.53	9.7	4.66	40.9	
P3041-04	FT011P-SD04-062024-00	10	1.15	8.83	9.98	4.49	37.8	
P3041-05	FT013P-SD03-062024-00	11	1.18	8.36	9.54	5.8	55.3	
P3041-06	AT029P-SD04-062024-00	12	1.17	8.53	9.7	7.74	77.0	
P3041-07	SS030P-SD02-062024-00	13	1.19	8.79	9.98	5.89	53.5	
P3041-09	AT029P-SD01-062024-00	14	1.17	8.47	9.64	6.86	67.2	
P3041-10	P3041-09MS	15	1.17	8.47	9.64	6.86	67.2	
P3041-11	P3041-09MSD	16	1.17	8.47	9.64	6.86	67.2	
P3041-12	DA053P-SD03-062024-00	17	1.17	8.45	9.62	6.99	68.9	
P3042-01	SS055P-SD04-062024-00	18	1.17	8.72	9.89	7.95	77.8	
P3042-02	SS057P-SD01-062024-00	19	1.17	8.57	9.74	7.18	70.1	
P3042-03	WP021P-SD01-062024-00	20	1.18	8.80	9.98	3.5	26.4	
P3042-04	SS057P-SD03-062024-00	21	1.17	8.80	9.97	7.29	69.5	
P3042-05	SS057P-SD02-062024-00	22	1.16	8.82	9.98	3.19	23.0	
P3042-06	SS057P-SD01-062024-00	23	1.17	8.43	9.6	6.62	64.7	
P3043-01	SS056P-SD05-062024-00	24	1.19	8.79	9.98	7.81	75.3	
P3043-02	SS056P-SD04-062024-00	25	1.18	8.79	9.97	7.14	67.8	
P3043-03	SS052P-SD03-062024-00	26	1.17	8.36	9.53	7.06	70.5	
P3043-04	SS056P-SD01-062024-00	27	1.17	8.79	9.96	8.58	84.3	
P3043-05	SS052P-SD02-062024-00	28	1.18	8.55	9.73	7.33	71.9	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/25/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 06/24/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/25/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB131387

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3045-01	6AB	29	1.18	8.47	9.65	9.44	97.5	
P3045-03	7AB	30	1.14	8.44	9.58	9.56	99.8	
P3045-05	8AB	31	1.15	8.38	9.53	9.4	98.4	
P3045-07	9AB	32	1.15	8.38	9.53	9.43	98.8	
P3045-09	10AB	33	1.18	8.72	9.9	9.85	99.4	
P3046-01	KMA963J-1-1	34	1.00	1.00	2.00	2.00	100.0	oilc
P3046-02	KMA963J-1-2	35	1.00	1.00	2.00	2.00	100.0	oilc
P3046-03	BC260098-1-1	36	1.00	1.00	2.00	2.00	100.0	oilc
P3046-04	BC260098-1-2	37	1.00	1.00	2.00	2.00	100.0	oilc
P3047-01	EO-01-062424	38	1.15	8.37	9.52	9.1	95.0	
P3047-02	EO-01-062424-E2	39	1.12	8.76	9.88	9.35	93.9	
P3048-01	Y2308-0064-1-1	40	1.00	1.00	2.00	2.00	100.0	oilc
P3048-02	Y2308-0064-1-2	41	1.00	1.00	2.00	2.00	100.0	oilc
P3048-03	Y2308-0068-1-1	42	1.00	1.00	2.00	2.00	100.0	oilc
P3048-04	Y2308-0068-1-2	43	1.00	1.00	2.00	2.00	100.0	oilc
P3048-05	BC271290-1-1	44	1.00	1.00	2.00	2.00	100.0	oilc
P3048-06	BC271290-1-2	45	1.00	1.00	2.00	2.00	100.0	oilc
P3048-07	BC288016-1-1	46	1.00	1.00	2.00	2.00	100.0	oilc
P3048-08	BC288016-1-2	47	1.00	1.00	2.00	2.00	100.0	oilc
P3048-09	JMA219R-1-1	48	1.00	1.00	2.00	2.00	100.0	oilc
P3048-10	JMA219R-1-2	49	1.00	1.00	2.00	2.00	100.0	oilc
P3048-11	JEC425W-1-1	50	1.00	1.00	2.00	2.00	100.0	oilc
P3048-12	JEC425W-1-2	51	1.00	1.00	2.00	2.00	100.0	oilc
P3049-01	SHOP-RAG-DRUM	52	1.00	1.00	2.00	2.00	100.0	debris
P3050-01	624-01	53	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-02	624-02	54	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-03	624-03	55	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-04	624-04	56	1.00	1.00	2.00	2.00	100.0	wipe sample

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/25/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 06/24/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/25/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB131387

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3050-05	624-05	57	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-06	624-06	58	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-07	624-07	59	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-08	624-08	60	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-09	624-09	61	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-10	624-10	62	1.00	1.00	2.00	2.00	100.0	wipe sample
P3050-11	624-11	63	1.00	1.00	2.00	2.00	100.0	wipe sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$