

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/28/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/27/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 06/28/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131428

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3098-01	AT014P-SD07-062524-00	1	1.18	8.66	9.84	7.94	78.1	
P3098-02	AT014P-SD04-062524-00	2	1.15	8.72	9.87	8.12	79.9	
P3098-03	SS092P-SD06-062524-00	3	1.19	8.79	9.98	9.00	88.9	
P3098-04	SS092P-SD07-062524-00	4	1.17	8.77	9.94	8.05	78.4	
P3098-05	SS091P-SD05-062524-00	5	1.15	8.54	9.69	8.37	84.5	
P3098-06	AT090P-SD06-062524-00	6	1.16	8.78	9.94	7.19	68.7	
P3098-07	AT016P-SD05-062524-00	7	1.17	8.80	9.97	7.5	71.9	
P3106-01	S-1	8	1.15	8.43	9.58	9.54	99.5	
P3106-02	S-2	9	1.18	8.40	9.58	9.54	99.5	
P3106-03	S-3	10	1.16	8.47	9.63	9.59	99.5	
P3106-04	FILTER-1	11	1.16	8.48	9.64	9.41	97.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$