

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 6/28/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 06/27/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:44
Out Date: 06/28/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131432

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3033-01	MA4CF0	1	1.18	8.68	9.86	8.61	85.6	
P3033-02	MA4CG2	2	1.15	8.38	9.53	9.11	95.0	
P3033-03	MA4CG2D	3	1.15	8.38	9.53	9.11	95.0	
P3033-04	MA4CG2S	4	1.15	8.38	9.53	9.11	95.0	
P3033-05	MA4CG3	5	1.15	8.52	9.67	9.36	96.4	
P3033-06	MA4B96	6	1.17	8.82	9.99	9.38	93.1	
P3033-07	MA4B97	7	1.19	8.78	9.97	8.73	85.9	
P3033-08	MA4CB8	8	1.13	8.81	9.94	7.78	75.5	
P3033-09	MA4CB9	9	1.15	8.80	9.95	8.5	83.5	
P3033-11	MA4D72	10	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3033-12	MA4CC0	11	1.16	8.67	9.83	9.3	93.9	
P3033-13	MA4CC4	12	1.17	8.81	9.98	9.22	91.4	
P3033-14	MA4CC9	13	1.17	8.82	9.99	8.12	78.8	
P3033-15	MA4CD0	14	1.17	8.81	9.98	8.6	84.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$