

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 7/2/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 07/01/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 07/02/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131467

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P2983-01	PCB061424RB02	1	1.00	1.00	2.00	2.00	100.0	caluk
P3023-02	S01-01C	9	1.18	8.28	9.46	8.99	94.3	
P3023-04	S01-02C	10	1.15	8.51	9.66	9.26	95.3	
P3023-06	S01-03C	11	1.18	8.56	9.74	8.83	89.4	
P3023-08	S02-01C	12	1.12	8.76	9.88	9.00	90.0	
P3023-10	S02-02C	13	1.18	8.56	9.74	8.73	88.2	
P3023-12	S02-03C	14	1.15	8.76	9.91	8.78	87.1	
P3023-14	S02-04C	15	1.19	8.67	9.86	9.06	90.8	
P3023-16	S02-05C	16	1.19	8.51	9.7	8.62	87.3	
P3023-18	S02-06C	17	1.16	8.48	9.64	8.66	88.4	
P3141-01	B-140-SB00	2	1.18	8.78	9.96	9.09	90.1	
P3141-02	B-140-SB01	3	1.19	8.65	9.84	8.41	83.5	
P3141-03	B-140-SB02	4	1.15	8.84	9.99	7.3	69.6	
P3141-04	B-129-SB01	5	1.18	8.76	9.94	8.7	85.8	
P3141-05	B-129-SB02	6	1.12	8.77	9.89	6.69	63.5	
P3142-01	VNJ-231	7	1.00	1.00	2.00	2.00	100.0	conc sample,100 & solids
P3142-03	327	8	1.00	1.00	2.00	2.00	100.0	conc sample,100 & solids

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$