

PERCENT SOLID

Supervisor: Iwona
Analyst: Niha
Date: 7/16/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 13:10
In Date: 07/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN #1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 07/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131627

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3216-01	MA4C00	1	1.14	8.45	9.59	7.67	77.3	
P3216-02	MA4C02	2	1.16	8.74	9.9	8.74	86.7	
P3216-03	MA4C08	3	1.12	8.62	9.74	8.8	89.1	
P3216-04	MA4C09	4	1.16	8.43	9.59	8.43	86.2	
P3216-05	MA4C12	5	1.16	8.62	9.78	8.91	89.9	
P3216-06	MA4D23	6	1.19	8.52	9.71	3.85	31.2	
P3216-07	MA4D23D	7	1.19	8.52	9.71	3.85	31.2	
P3216-08	MA4D23S	8	1.19	8.52	9.71	3.85	31.2	
P3216-09	MA4D24	9	1.15	8.38	9.53	4.28	37.4	
P3216-10	MA4D25	10	1.18	8.63	9.81	4.45	37.9	
P3216-11	MA4BD3	11	1.14	8.47	9.61	8.64	88.5	
P3216-12	MA4BD4	12	1.16	8.71	9.87	8.82	87.9	
P3216-13	MA4D26	13	1.17	8.56	9.73	9.26	94.5	
P3216-14	MA4D31	14	1.17	8.38	9.55	8.77	90.7	
P3216-15	MA4D43	15	1.19	8.55	9.74	7.65	75.6	
P3216-16	MA4D44	16	1.17	8.74	9.91	8.32	81.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$