

PERCENT SOLID

Supervisor: Iwona
Analyst: Niha
Date: 7/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 07/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 07/18/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131649

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3288-01	1	1	1.00	1.00	2.00	2.00	100.0	caluk
P3288-02	2	2	1.00	1.00	2.00	2.00	100.0	caluk
P3288-03	3	3	1.00	1.00	2.00	2.00	100.0	caluk
P3288-04	4	4	1.00	1.00	2.00	2.00	100.0	caluk
P3288-05	5	5	1.00	1.00	2.00	2.00	100.0	caluk
P3291-01	PIPE-9	6	1.00	1.00	2.00	2.00	100.0	wipe sample
P3293-01	716-A	7	1.00	1.00	2.00	2.00	100.0	wipe sample
P3293-02	716-B	8	1.00	1.00	2.00	2.00	100.0	wipe sample
P3293-03	716-C	9	1.00	1.00	2.00	2.00	100.0	wipe sample
P3294-01	HR-1-071724	10	1.18	8.59	9.77	9.54	97.3	
P3294-02	HR-1-071724-E2	11	1.2	8.54	9.74	9.62	98.6	
P3294-03	HR-3-071724	12	1.19	8.63	9.82	9.48	96.1	
P3294-04	HR-3-071724-E2	13	1.18	8.67	9.85	9.52	96.2	
P3294-05	HR-4-071724	14	1.17	8.55	9.72	8.99	91.5	
P3294-06	HR-4-071724-E2	15	1.17	8.78	9.95	9.03	89.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$