

PERCENT SOLID

Supervisor: Iwona
Analyst: Niha
Date: 7/22/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 07/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 07/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131682

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3283-01	A4CM2	1	1.18	8.60	9.78	9.38	95.3	
P3283-02	A4DJ0	2	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3283-03	A4DQ9	3	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3283-04	A4DH7	4	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3283-05	VHBLK001	5	1.00	1.00	2.00	2.00	100.0	vhblk
P3283-06	A4CM2MS	6	1.18	8.60	9.78	9.38	95.3	
P3283-07	A4CM2MSD	7	1.18	8.60	9.78	9.38	95.3	
P3285-01	MA4CM2	8	1.18	8.60	9.78	9.38	95.3	
P3285-02	MA4D97	9	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3285-03	MA4CM2D	10	1.18	8.60	9.78	9.38	95.3	
P3285-04	MA4CM2S	11	1.18	8.60	9.78	9.38	95.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$