

PERCENT SOLID

Supervisor: Iwona
Analyst: Niha
Date: 7/22/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 14:30
In Date: 07/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 07/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131683

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3265-01	A4DH6	1	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3265-02	A4DL0	2	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3265-03	A4DQ8	3	1.00	1.00	2.00	2.00	100.0	P.T. SAMPLE
P3265-04	VHBLK001	4	1.00	1.00	2.00	2.00	100.0	vhblk
P3265-05	A4D22	5	1.18	8.71	9.89	6.68	63.1	
P3265-06	A4D39	6	1.18	8.63	9.81	8.17	81.0	
P3265-07	A4D40	7	1.18	8.58	9.76	8.3	83.0	
P3265-08	A4D41	8	1.18	8.80	9.98	6.42	59.5	
P3265-09	A4D42	9	1.18	8.34	9.52	7.53	76.1	
P3265-10	A4DA8	10	1.00	1.00	2.00	2.00	100.0	FB
P3265-11	A4D42MS	11	1.18	8.34	9.52	7.53	76.1	
P3265-12	A4D42MSD	12	1.18	8.34	9.52	7.53	76.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$