

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131878

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3474-01	WC-A1-COMP	1	1.18	8.64	9.82	8.8	88.2	
P3474-03	WC-A1-GRAB	2	1.15	8.38	9.53	8.78	91.1	
P3474-04	WC-B1-COMP	3	1.16	8.59	9.75	8.56	86.1	
P3474-06	WC-B1-GRAB	4	1.18	8.43	9.61	8.18	83.0	
P3474-07	WC-A2-COMP	5	1.19	8.65	9.84	8.25	81.6	
P3474-09	WC-A2-GRAB	6	1.16	8.69	9.85	8.12	80.1	
P3474-10	WC-B2-COMP	7	1.16	8.67	9.83	8.25	81.8	
P3474-12	WC-B2-GRAB	8	1.16	8.69	9.85	8.09	79.7	
P3475-01	S-866-N-SO-0-0.5-080524	9	1.14	8.83	9.97	8.49	83.2	
P3475-02	S-866-N-SO-1.0-1.5-080524	10	1.15	8.66	9.81	8.57	85.7	
P3475-03	P3475-02MS	11	1.15	8.66	9.81	8.57	85.7	
P3475-04	P3475-02MSD	12	1.15	8.66	9.81	8.57	85.7	
P3475-05	S-858-N-SO-1.0-1.5-080524	13	1.12	8.74	9.86	7.56	73.7	sludge sample
P3475-06	P3475-05MS	14	1.12	8.74	9.86	7.56	73.7	sludge sample
P3475-07	P3475-05MSD	15	1.12	8.74	9.86	7.56	73.7	sludge sample
P3476-02	MUL04	16	1.15	8.50	9.65	5.16	47.2	
P3476-04	MUL05	17	1.16	8.81	9.97	4.97	43.2	
P3476-06	MUL06	18	1.15	8.75	9.9	4.98	43.8	
P3476-08	MUL07	19	1.14	8.60	9.74	5.07	45.7	
P3478-04	TS04	20	1.18	8.57	9.75	8.17	81.6	
P3478-06	TS05	21	1.16	8.61	9.77	8.59	86.3	
P3478-08	TS06	22	1.15	8.83	9.98	8.45	82.7	
P3478-10	TS07	23	1.14	8.83	9.97	8.52	83.6	
P3479-02	TS08	24	1.14	8.42	9.56	7.99	81.4	
P3479-04	TS09	25	1.14	8.71	9.85	7.98	78.5	
P3479-06	TS10	26	1.14	8.58	9.72	8.15	81.7	
P3479-08	TS11	27	1.15	8.49	9.64	8.79	90.0	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131878

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3483-01	A	28	1.12	8.46	9.58	9.11	94.4	
P3484-01	B	29	1.18	8.72	9.9	9.22	92.2	
P3485-01	C	30	1.11	8.88	9.99	8.64	84.8	
P3487-01	0726-A	31	1.00	1.00	2.00	2.00	100.0	oil sample
P3487-02	0726-B	32	1.00	1.00	2.00	2.00	100.0	oil sample
P3487-03	0726-C	33	1.00	1.00	2.00	2.00	100.0	oil sample
P3487-05	0723-B	34	1.00	1.00	2.00	2.00	100.0	oil sample
P3488-01	HD-01-080624	35	1.14	8.58	9.72	8.88	90.2	
P3488-02	HD-01-080624-E2	36	1.15	8.80	9.95	8.96	88.8	
P3488-03	HD-02-080624	37	1.14	8.40	9.54	8.76	90.7	
P3488-04	HD-02-080624-E2	38	1.15	8.81	9.96	9.34	93.0	
P3489-01	WASTE	39	1.15	8.62	9.77	5.9	55.1	
P3489-02	VOC	40	1.19	8.52	9.71	5.72	53.2	
P3489-03	1	41	1.16	8.81	9.97	6.32	58.6	
P3489-04	2	42	1.16	8.81	9.97	5.57	50.1	
P3489-05	3	43	1.15	8.38	9.53	5.38	50.5	
P3489-06	4	44	1.19	8.42	9.61	6.12	58.6	
P3489-07	5	45	1.15	8.80	9.95	6.02	55.3	
P3491-01	72-11936	46	1.18	8.54	9.72	9.47	97.1	
P3491-02	SOM-24-00087	47	1.00	1.00	2.00	2.00	100.0	debris
P3491-03	COMP-1	48	1.00	1.00	2.00	2.00	100.0	debris
P3491-04	COMP-2	49	1.00	1.00	2.00	2.00	100.0	debris
P3491-05	SOM-24-00094	50	1.00	1.00	2.00	2.00	100.0	oily-debris
P3492-01	Y2309-0098-1-1	51	1.00	1.00	2.00	2.00	100.0	oilc
P3492-02	Y2309-0098-1-2	52	1.00	1.00	2.00	2.00	100.0	oilc
P3493-01	628-A	53	1.17	8.41	9.58	9.54	99.5	
P3493-02	628-B	54	1.16	8.42	9.58	9.51	99.2	
P3494-01	ORA-1899	55	1.00	1.00	2.00	2.00	100.0	wipe sample



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 08/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131878

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3494-02	ORA-1900	56	1.00	1.00	2.00	2.00	100.0	wipe sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$