



PERCENT SOLID

Supervisor: Sohil
Analyst: jignesh
Date: 8/8/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 08/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:50
Out Date: 08/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB131895

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3506-01	RT4534	1	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P3507-01	EO-03-080724	2	1.16	8.48	9.64	8.95	91.9	
P3507-02	EO-03-080724-E2	3	1.16	8.65	9.81	8.97	90.3	
P3508-01	BC260544	4	1.00	1.00	2.00	2.00	100.0	pilc
P3508-02	BC260544	5	1.00	1.00	2.00	2.00	100.0	pilc
P3508-03	KMH2650	6	1.00	1.00	2.00	2.00	100.0	pilc
P3508-04	KMH2650	7	1.00	1.00	2.00	2.00	100.0	pilc
P3508-05	BC288026	8	1.00	1.00	2.00	2.00	100.0	pilc
P3508-06	BC288026	9	1.00	1.00	2.00	2.00	100.0	pilc
P3508-07	TERMINATOR	10	1.00	1.00	2.00	2.00	100.0	pilc
P3509-01	OR-02-080724	11	1.15	8.76	9.91	8.32	81.8	
P3509-02	OR-02-080724-E2	12	1.13	8.59	9.72	7.97	79.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$