

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 8/16/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 14:05
In Date: 08/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 08/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132030

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3571-01	E1JD8	1	1.14	8.58	9.72	7.95	79.4	
P3571-02	E1JD9	2	1.15	8.55	9.7	7.87	78.6	
P3571-03	E1JE0	3	1.18	8.48	9.66	8.03	80.8	
P3571-04	E1JE0D	4	1.18	8.48	9.66	8.03	80.8	
P3571-05	E1JE0S	5	1.18	8.48	9.66	8.03	80.8	
P3571-06	E1JM2	6	1.13	8.55	9.68	8.01	80.5	
P3571-07	E1JM3	7	1.18	8.72	9.9	8.00	78.2	
P3571-08	E1JM4	8	1.15	8.84	9.99	8.13	79.0	
P3571-09	E1JM5	9	1.12	8.70	9.82	8.23	81.7	
P3571-10	E1JM6	10	1.15	8.84	9.99	7.98	77.3	
P3571-11	E1JM7	11	1.18	8.66	9.84	7.82	76.7	
P3571-12	E1JM8	12	1.15	8.82	9.97	8.12	79.0	
P3571-13	E1JM9	13	1.15	8.83	9.98	8.62	84.6	
P3571-14	E1JN0	14	1.14	8.81	9.95	8.32	81.5	
P3571-15	E1JN1	15	1.19	8.57	9.76	7.91	78.4	
P3571-16	E1JN2	16	1.17	8.60	9.77	8.32	83.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$