

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/4/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 09/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132252

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3673-09	DCXX9	1	1.15	8.53	9.68	7.9	79.1	
P3816-01	SS050P-SD23-082824-00	2	1.15	8.81	9.96	7.75	74.9	
P3816-02	SS050P-SD24-082824-00	3	1.15	8.84	9.99	7.48	71.6	
P3816-03	SS050P-SD25-082824-00	4	1.15	8.82	9.97	7.75	74.8	
P3816-04	SS050P-SD26-082824-00	5	1.12	8.47	9.59	7.09	70.5	
P3816-05	SS050P-SD20-082824-00	6	1.15	8.82	9.97	9.00	89.0	
P3816-06	SS050P-SD22-082924-00	7	1.18	8.50	9.68	5.46	50.4	
P3816-07	SS050P-SD21-082924-00	8	1.15	8.55	9.7	7.24	71.2	
P3816-08	SS055P-SD20-082924-00	9	1.14	8.82	9.96	5.74	52.2	
P3816-09	SS054P-SD30-082924-00	10	1.14	8.60	9.74	8.63	87.1	
P3816-10	SS054P-SD35-082924-00	11	1.12	8.77	9.89	8.33	82.2	
P3816-11	SS054P-SD36-082924-00	12	1.13	8.54	9.67	2.35	14.3	
P3816-12	SS055P-SD06-082924-00	13	1.15	8.67	9.82	5.49	50.1	
P3817-01	AT024P-SD40-082724-00	14	1.18	8.81	9.99	8.58	84.0	
P3817-02	AT024P-SD29-082724-00	15	1.18	8.59	9.77	6.24	58.9	
P3817-03	AT024P-SD30-082724-00	16	1.14	8.67	9.81	4.14	34.6	
P3817-04	AT024P-SD38-082724-00	17	1.12	8.86	9.98	7.67	73.9	
P3817-05	AT024P-SD39-082724-00	18	1.12	8.65	9.77	7.66	75.6	
P3817-06	AT024P-SD42-082724-00	19	1.12	8.72	9.84	6.95	66.9	
P3817-07	AT024P-SD22-082724-00	20	1.15	8.83	9.98	6.52	60.8	
P3817-08	AT024P-SD35-082824-00	21	1.15	8.72	9.87	3.1	22.4	
P3817-09	AT024P-SD36-082824-00	22	1.15	8.82	9.97	6.94	65.6	
P3818-01	AT014P-SD09-082924-00	23	1.15	8.82	9.97	8.11	78.9	
P3818-02	AT014P-SD02-082924-00	24	1.18	8.46	9.64	7.73	77.4	
P3818-03	AT090P-SD01-083024-00	25	1.15	8.83	9.98	8.13	79.0	
P3818-04	AT090P-SD02-083024-00	26	1.16	8.80	9.96	7.8	75.5	
P3818-05	SS092P-SD01-083024-00	27	1.15	8.36	9.51	7.83	79.9	
P3818-06	SS091P-SD02-083024-00	28	1.14	8.81	9.95	8.78	86.7	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/4/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 09/03/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132252

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3818-07	SS091PP-SD01-083024-00	29	1.15	8.82	9.97	7.4	70.9	
P3818-08	AT016P-SD01-083024-00	30	1.17	8.60	9.77	6.76	65.0	
P3818-09	AT016P-SD02-083024-00	31	1.15	8.80	9.95	2.85	19.3	
P3818-10	SS091P-SD03-083024-00	32	1.15	8.83	9.98	8.1	78.7	
P3819-01	OCB-CONCRETE	33	1.00	1.00	2.00	2.00	100.0	conc sample,100 % solids
P3819-02	OCB-CONCRETE	34	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P3820-01	STONE-1	35	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P3820-02	STONE-1	36	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P3822-01	WE-1	37	1.00	1.00	2.00	2.00	100.0	wipe sample
P3822-02	WE-2	38	1.00	1.00	2.00	2.00	100.0	wipe sample
P3825-01	DN-W-01 (B-1)	39	1.13	8.75	9.88	8.51	84.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$