

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/6/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 09/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 09/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132290

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3849-01	MP-1	1	1.00	1.00	2.00	2.00	100.0	wipe sample
P3849-02	MP-2	2	1.00	1.00	2.00	2.00	100.0	wipe sample
P3849-03	MP-3	3	1.00	1.00	2.00	2.00	100.0	wipe sample
P3849-04	MP-4	4	1.00	1.00	2.00	2.00	100.0	wipe sample
P3849-05	MP-5	5	1.00	1.00	2.00	2.00	100.0	wipe sample
P3851-01	BC122669-1-1	6	1.00	1.00	2.00	2.00	100.0	oilc
P3851-02	BC122669-1-2	7	1.00	1.00	2.00	2.00	100.0	oilc
P3851-03	BC271332-1-1	8	1.00	1.00	2.00	2.00	100.0	oilc
P3851-04	BC271332-1-2	9	1.00	1.00	2.00	2.00	100.0	oilc
P3851-05	BC271332-2-1	10	1.00	1.00	2.00	2.00	100.0	oilc
P3851-06	BC271332-2-2	11	1.00	1.00	2.00	2.00	100.0	oilc
P3851-07	BC238813-1-1	12	1.00	1.00	2.00	2.00	100.0	oilc
P3851-08	BC238813-1-2	13	1.00	1.00	2.00	2.00	100.0	oilc
P3851-09	BC238813-2-1	14	1.00	1.00	2.00	2.00	100.0	oilc
P3851-10	BC238813-2-2	15	1.00	1.00	2.00	2.00	100.0	oilc
P3851-11	BC271246-1-1	16	1.00	1.00	2.00	2.00	100.0	oilc
P3851-12	BC271246-1-2	17	1.00	1.00	2.00	2.00	100.0	oilc
P3851-13	BC274767-1-1	18	1.00	1.00	2.00	2.00	100.0	oilc
P3851-14	BC274767-1-2	19	1.00	1.00	2.00	2.00	100.0	oilc
P3852-01	D3656	20	1.00	1.00	2.00	2.00	100.0	debris
P3852-02	VNJ223	21	1.15	8.64	9.79	9.29	94.2	
P3852-03	VNJ223-E2	22	1.15	8.75	9.9	9.39	94.2	
P3852-04	ETGI356	23	1.15	8.57	9.72	8.00	79.9	
P3852-05	D356-E2	24	1.17	8.58	9.75	8.93	90.4	
P3854-01	40295	25	1.00	1.00	2.00	2.00	100.0	wipe sample
P3854-02	40296	26	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-01	PIPE-1	27	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-02	PIPE-2	28	1.00	1.00	2.00	2.00	100.0	wipe sample



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/6/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 09/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 09/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132290

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3855-03	PIPE-3	29	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-04	PIPE-4	30	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-05	PIPE-5	31	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-06	PIPE-6	32	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-07	PIPE-7	33	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-08	PIPE-8	34	1.00	1.00	2.00	2.00	100.0	wipe sample
P3855-09	PIPE-9	35	1.00	1.00	2.00	2.00	100.0	wipe sample
P3857-01	WC-C1-COMP	36	1.14	8.78	9.92	8.31	81.7	
P3857-03	WC-C1-GRAB	37	1.15	8.40	9.55	8.07	82.4	
P3857-04	WC-C2-COMP	38	1.18	8.48	9.66	8.07	81.3	
P3857-06	WC-C2-GRAB	39	1.15	8.83	9.98	8.58	84.1	
P3857-07	WC-D1-COMP	40	1.15	8.50	9.65	8.33	84.5	
P3857-09	WC-D1-GRAB	41	1.15	8.82	9.97	8.73	85.9	
P3857-10	WC-D2-COMP	42	1.15	8.68	9.83	8.31	82.5	
P3857-12	WC-D2-GRAB	43	1.16	8.43	9.59	8.17	83.2	
P3861-01	DN-B-37	44	1.18	8.57	9.75	9.36	95.4	
P3866-01	MR-CAM-04	45	1.00	1.00	2.00	2.00	100.0	STONES

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$