



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/9/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:28
Out Date: 09/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132307

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P3856-01	WASTE	1	1.12	8.77	9.89	9.52	95.8	
P3856-02	VOC	2	1.12	8.70	9.82	9.46	95.9	
P3856-03	1	3	1.15	8.40	9.55	9.2	95.8	
P3856-04	2	4	1.18	8.75	9.93	9.45	94.5	
P3856-05	3	5	1.18	8.65	9.83	9.4	95.0	
P3856-06	4	6	1.13	8.65	9.78	9.45	96.2	
P3856-07	5	7	1.12	8.60	9.72	9.38	96.0	
P3871-01	B1 (0-6)	18	1.18	8.47	9.65	8.41	85.4	
P3876-05	SVOC-GPC-BLANK	8	1.00	1.00	2.00	2.00	100.0	
P3876-06	PEST-GPC-BLANK	9	1.00	1.00	2.00	2.00	100.0	
P3876-07	PEST-GPC-BLANK-SPIKE	10	1.00	1.00	2.00	2.00	100.0	
P3876-08	PCB-GPC-BLANK	11	1.00	1.00	2.00	2.00	100.0	
P3876-09	PCB-GPC-BLANK-SPIKE	12	1.00	1.00	2.00	2.00	100.0	
P3876-10	SVOC-GPC2-BLANK	13	1.00	1.00	2.00	2.00	100.0	
P3876-11	PEST-GPC2-BLANK	14	1.00	1.00	2.00	2.00	100.0	
P3876-12	PEST-GPC2-BLANK-SPIKE	15	1.00	1.00	2.00	2.00	100.0	
P3876-13	PCB-GPC2-BLANK	16	1.00	1.00	2.00	2.00	100.0	
P3876-14	PCB-GPC2-BLANK-SPIKE	17	1.00	1.00	2.00	2.00	100.0	
P3880-01	SHEEPSHEAD-SOIL-PILE	19	1.15	8.42	9.57	8.83	91.2	
P3880-03	SHEEPSHEAD-SOIL-PILE	20	1.15	8.75	9.9	9.09	90.7	
P3881-01	VNJ-217	21	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P3887-01	60201	22	1.00	1.00	2.00	2.00	100.0	wipe sample
P3887-02	60202	23	1.00	1.00	2.00	2.00	100.0	wipe sample
P3888-01	DN-B-36	24	1.14	8.78	9.92	9.22	92.0	
P3888-02	DN-W-02	25	1.12	8.62	9.74	8.63	87.1	
P3892-01	ARS20-005	26	1.15	8.84	9.99	8.72	85.6	
P3893-01	KMA935G-1-1	27	1.00	1.00	2.00	2.00	100.0	pile
P3893-03	KMA935G-1-2	28	1.00	1.00	2.00	2.00	100.0	pile



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/9/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:28
Out Date: 09/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132307

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P3894-01	SHOP-DRUM	29	1.18	8.71	9.89	5.96	54.9	
P3895-01	FDA6470-1-1	30	1.00	1.00	2.00	2.00	100.0	oilc
P3895-02	FDA6470-1-2	31	1.00	1.00	2.00	2.00	100.0	oilc
P3896-02	MOO-23-00293	32	1.13	8.67	9.8	9.33	94.6	
P3897-01	1	33	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-02	2	34	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-03	3	35	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-04	4	36	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-05	5	37	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-06	6	38	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-07	7	39	1.00	1.00	2.00	2.00	100.0	wipe sample
P3897-08	8	40	1.00	1.00	2.00	2.00	100.0	wipe sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$