



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 09/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132351

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P3923-01	OILY-STONE	1	1.00	1.00	2.00	2.00	100.0	oily stone
P3925-01	NB-07091024	2	1.18	8.34	9.52	9.12	95.2	
P3925-02	NB-07091024-E2	3	1.18	8.34	9.52	9.03	94.1	
P3928-01	091024-A	4	1.00	1.00	2.00	2.00	100.0	wipe sample
P3928-02	091024-B	5	1.00	1.00	2.00	2.00	100.0	wipe sample
P3928-03	091024-C	6	1.00	1.00	2.00	2.00	100.0	wipe sample
P3928-04	091024-D	7	1.00	1.00	2.00	2.00	100.0	wipe sample
P3929-01	BIN-2	8	1.15	8.40	9.55	8.87	91.9	
P3929-02	BIN-2	9	1.18	8.43	9.61	8.7	89.2	
P3929-03	72-11930	10	1.13	8.74	9.87	8.65	86.0	
P3929-04	72-11930	11	1.13	8.70	9.83	8.97	90.1	
P3930-01	DN-B-33	12	1.15	8.43	9.58	8.65	89.0	
P3930-02	DN-B-34	13	1.18	8.48	9.66	9.18	94.3	
P3932-01	FM1	14	1.17	8.80	9.97	9.33	92.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$