



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/13/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 15:25
In Date: 09/12/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 09/13/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132406

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P3974-01	EIKS8	1	1.18	8.66	9.84	8.9	89.1	
P3974-02	EIKT0	2	1.15	8.43	9.58	8.02	81.5	
P3974-03	EIKT2	3	1.15	8.83	9.98	8.76	86.2	
P3974-04	EIKT4	4	1.15	8.38	9.53	8.22	84.4	
P3974-05	EIKT6	5	1.17	8.60	9.77	8.42	84.3	
P3974-06	EIKU3	6	1.18	8.52	9.7	8.59	87.0	
P3974-07	EIKU5	7	1.15	8.82	9.97	8.82	87.0	
P3974-08	EIKU7	8	1.18	8.62	9.8	9.39	95.2	
P3974-09	EIKU9	9	1.15	8.76	9.91	9.54	95.8	
P3974-10	EIKV1	10	1.14	8.40	9.54	9.22	96.2	
P3974-11	EIKV4	11	1.15	8.40	9.55	8.38	86.1	
P3974-12	EIKV6	12	1.13	8.80	9.93	9.48	94.9	
P3974-13	EIKV8	13	1.17	8.59	9.76	9.42	96.0	
P3974-14	EIKW0	14	1.15	8.48	9.63	9.26	95.6	
P3974-15	EIKW2	15	1.15	8.44	9.59	9.25	96.0	
P3974-16	EIKW4	16	1.18	8.43	9.61	9.27	96.0	
P3974-17	EIKW4D	17	1.18	8.43	9.61	9.27	96.0	
P3974-18	EIKW4S	18	1.18	8.43	9.61	9.27	96.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$