



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/16/2024

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 09/13/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 09/14/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132440

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4006-01	MH2A01	1	1.15	8.83	9.98	7.18	68.3	
P4006-02	MH2A02	2	1.18	8.63	9.81	8.12	80.4	
P4006-03	MH2A03	3	1.16	8.50	9.66	7.02	68.9	
P4006-04	MH2A04	4	1.15	8.76	9.91	8.14	79.8	
P4006-05	MH2A05	5	1.14	8.81	9.95	7.46	71.7	
P4006-06	MH2A06	6	1.15	8.81	9.96	8.32	81.4	
P4006-07	MH2A07	7	1.15	8.80	9.95	8.25	80.7	
P4006-08	MH2A08	8	1.17	8.60	9.77	7.16	69.7	
P4006-09	MH2A09	9	1.18	8.50	9.68	7.08	69.4	
P4006-10	MH2A10	10	1.15	8.58	9.73	8.17	81.8	
P4006-11	MH2A11	11	1.18	8.49	9.67	7.57	75.3	
P4006-12	MH2A12	12	1.18	8.61	9.79	7.79	76.8	
P4006-13	MH2A13	13	1.15	8.82	9.97	8.78	86.5	
P4006-14	MH2A14	14	1.11	8.66	9.77	7.5	73.8	
P4006-15	MH2A15	15	1.15	8.80	9.95	8.67	85.5	
P4006-16	MH2A16	16	1.16	8.40	9.56	8.18	83.6	
P4006-17	MH2A17	17	1.15	8.47	9.62	7.13	70.6	
P4006-18	MH2A18	18	1.14	8.82	9.96	7.9	76.6	
P4006-19	MH2A18D	19	1.14	8.82	9.96	7.9	76.6	
P4006-20	MH2A18S	20	1.14	8.82	9.96	7.9	76.6	
P4006-21	MH2A19	21	1.14	8.80	9.94	7.71	74.7	
P4008-01	MH2A20	22	1.18	8.72	9.9	7.29	70.1	
P4008-02	MH2A20D	23	1.18	8.72	9.9	7.29	70.1	
P4008-03	MH2A20S	24	1.18	8.72	9.9	7.29	70.1	
P4008-04	MH2A21	25	1.13	8.80	9.93	7.65	74.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$