



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 09/18/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132467

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4085-01	COMP-1	1	1.12	8.65	9.77	9.02	91.3	
P4085-03	COMP-2	2	1.18	8.42	9.6	8.77	90.1	
P4086-01	SOIL	3	1.18	8.56	9.74	8.77	88.7	
P4086-02	SOIL-E2	4	1.15	8.59	9.74	8.56	86.3	
P4087-01	NB-614-COMP-01	15	1.18	8.42	9.6	8.65	88.7	
P4087-03	NB-614-VOC-01	16	1.18	8.46	9.64	8.57	87.4	
P4087-04	NB-614-01	17	1.12	8.76	9.88	8.74	87.0	
P4087-05	NB-614-02	18	1.16	8.67	9.83	8.69	86.9	
P4087-06	NB-614-03	19	1.12	8.60	9.72	8.31	83.6	
P4087-07	NB-614-COMP-02	20	1.15	8.41	9.56	8.35	85.6	
P4087-09	NB-614-VOC-02	21	1.12	8.72	9.84	8.81	88.2	
P4087-10	NB-614-04	22	1.13	8.68	9.81	8.34	83.1	
P4087-11	NB-614-05	23	1.12	8.45	9.57	8.14	83.1	
P4087-12	NB-614-06	24	1.16	8.77	9.93	8.9	88.3	
P4087-13	NB-614-COMP-03	25	1.18	8.48	9.66	8.43	85.5	
P4087-15	NB-614-VOC-03	26	1.14	8.40	9.54	7.77	78.9	
P4087-16	NB-614-07	27	1.12	8.75	9.87	8.61	85.6	
P4087-17	NB-614-08	28	1.15	8.84	9.99	9.04	89.3	
P4087-18	NB-614-09	29	1.13	8.77	9.9	8.49	83.9	
P4087-19	NB-614-COMP-04	30	1.15	8.81	9.96	8.69	85.6	
P4087-21	NB-614-VOC-04	31	1.12	8.70	9.82	8.34	83.0	
P4087-22	NB-614-10	32	1.18	8.59	9.77	8.42	84.3	
P4087-23	NB-614-11	33	1.15	8.83	9.98	8.36	81.7	
P4087-24	NB-614-12	34	1.15	8.82	9.97	8.75	86.2	
P4087-25	NB-614-COMP-05	35	1.18	8.76	9.94	8.3	81.3	
P4087-27	NB-614-VOC-05	36	1.13	8.84	9.97	7.02	66.6	
P4087-28	NB-614-13	37	1.18	8.74	9.92	7.29	69.9	
P4087-29	NB-614-14	38	1.18	8.43	9.61	7.85	79.1	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 09/18/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132467

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4087-30	NB-614-15	39	1.13	8.63	9.76	8.08	80.5	
P4087-31	NB-614-COMP-06	40	1.14	8.81	9.95	8.27	80.9	
P4087-33	NB-614-VOC-06	41	1.13	8.70	9.83	8.82	88.4	
P4087-34	NB-614-16	42	1.17	8.60	9.77	8.88	89.7	
P4087-35	NB-614-17	43	1.18	8.43	9.61	8.02	81.1	
P4087-36	NB-614-18	44	1.13	8.75	9.88	8.04	79.0	
P4088-01	NB-614-COMP-07	45	1.15	8.37	9.52	8.92	92.8	
P4088-03	NB-614-VOC-07	46	1.15	8.84	9.99	9.41	93.4	
P4088-04	NB-614-19	47	1.17	8.81	9.98	9.08	89.8	
P4088-05	NB-614-20	48	1.13	8.65	9.78	9.03	91.3	
P4088-06	NB-614-21	49	1.12	8.73	9.85	9.08	91.2	
P4088-07	NB-614-COMP-08	50	1.15	8.81	9.96	8.85	87.4	
P4088-09	NB-614-VOC-08	51	1.14	8.85	9.99	9.44	93.8	
P4088-10	NB-614-22	52	1.13	8.66	9.79	8.69	87.3	
P4088-11	NB-614-23	53	1.12	8.70	9.82	8.37	83.3	
P4088-12	NB-614-24	54	1.16	8.50	9.66	7.6	75.8	
P4088-13	NB-614-COMP-09	55	1.15	8.37	9.52	7.83	79.8	
P4088-15	NB-614-VOC-09	56	1.13	8.69	9.82	7.76	76.3	
P4088-16	NB-614-25	57	1.14	8.85	9.99	7.9	76.4	
P4088-17	NB-614-26	58	1.12	8.64	9.76	7.75	76.7	
P4088-18	NB-614-27	59	1.12	8.49	9.61	7.11	70.6	
P4089-01	365-351-294	11	1.15	8.84	9.99	8.57	83.9	
P4089-03	365-351-294-E2	12	1.15	8.62	9.77	7.82	77.4	
P4090-01	TR-04-091724	13	1.15	8.62	9.77	9.31	94.7	
P4090-02	TR-04-091724-E2	14	1.16	8.48	9.64	9.35	96.6	
P4097-01	RT2299	5	1.15	8.60	9.75	8.23	82.3	
P4097-02	RT2299	6	1.15	8.60	9.75	8.23	82.3	
P4098-01	HR-01-091724	7	1.18	8.58	9.76	9.24	93.9	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 9/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 09/18/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132467

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4098-02	HR-01-091724-E2	8	1.18	8.33	9.51	9.07	94.7	
P4098-03	HR-04-091724	9	1.15	8.40	9.55	9.2	95.8	
P4098-04	HR-04-091724E-E2	10	1.15	8.54	9.69	9.44	97.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$