

PERCENT SOLID

Supervisor: sohil
Analyst: jignesh
Date: 9/20/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132510

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4024-01	DDCD5	1	1.12	8.67	9.79	8.58	86.0	
P4024-02	DDCD5MS	2	1.12	8.67	9.79	8.58	86.0	
P4024-03	DDCD5MSD	3	1.12	8.67	9.79	8.58	86.0	
P4024-04	VHBLK001	4	1.00	1.00	2.00	2.00	100.0	vhblk
P4024-05	DDC36	5	1.16	8.40	9.56	7.77	78.7	
P4024-06	DDCL2	6	1.15	8.66	9.81	8.33	82.9	
P4103-01	WC-22A	7	1.18	8.49	9.67	9.05	92.7	
P4103-04	WC-14-5-7.5A	8	1.15	8.76	9.91	9.02	89.8	
P4103-07	WC-14-5-7.5B	9	1.14	8.83	9.97	9.1	90.1	
P4103-10	WC-14-7.5-10A	10	1.18	8.75	9.93	8.66	85.5	
P4103-13	WC-14-7.5-10B	11	1.14	8.42	9.56	8.56	88.1	
P4103-16	WC-23	12	1.15	8.73	9.88	8.66	86.0	
P4104-01	WC-14-5-7.5A-EPH-1	13	1.19	8.72	9.91	9.1	90.7	
P4104-02	WC-14-5-7.5A-EPH-2	14	1.17	8.61	9.78	8.95	90.4	
P4104-03	WC-14-5-7.5A-EPH-3	15	1.15	8.40	9.55	8.84	91.5	
P4104-04	WC-14-5-7.5A-EPH-4	16	1.12	8.84	9.96	8.48	83.3	
P4104-06	WC-14-5-7.5B-EPH-1	17	1.18	8.77	9.95	9.05	89.7	
P4104-07	WC-14-5-7.5B-EPH-2	18	1.15	8.62	9.77	8.87	89.6	
P4104-08	WC-14-5-7.5B-EPH-3	19	1.13	8.70	9.83	8.85	88.7	
P4104-09	WC-14-5-7.5B-EPH-4	20	1.15	8.82	9.97	8.91	88.0	
P4104-10	WC-14-5-7.5B-EPH-5	21	1.15	8.83	9.98	9.00	88.9	
P4104-11	WC-14-7.5-10A-EPH-1	22	1.15	8.80	9.95	9.19	91.4	
P4104-12	WC-14-7.5-10A-EPH-2	23	1.15	8.76	9.91	8.66	85.7	
P4104-13	WC-14-7.5-10A-EPH-3	24	1.15	8.81	9.96	9.12	90.5	
P4104-14	WC-14-7.5-10A-EPH-4	25	1.19	8.65	9.84	8.8	88.0	
P4104-15	WC-14-7.5-10A-EPH-5	26	1.15	8.41	9.56	8.7	89.8	
P4104-16	WC-14-7.5-10B-EPH-1	27	1.18	8.64	9.82	8.66	86.6	
P4104-17	WC-14-7.5-10B-EPH-2	28	1.15	8.83	9.98	8.75	86.1	

PERCENT SOLID

Supervisor: sohil
Analyst: jignesh
Date: 9/20/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132510

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4104-18	WC-14-7.5-10B-EPH-3	29	1.15	8.82	9.97	8.93	88.2	
P4104-19	WC-14-7.5-10B-EPH-4	30	1.19	8.57	9.76	8.76	88.3	
P4104-20	WC-14-7.5-10B-EPH-5	31	1.16	8.61	9.77	8.81	88.9	
P4104-21	WC-22A-EPH-1	32	1.15	8.81	9.96	9.44	94.1	
P4104-22	WC-22A-EPH-2	33	1.18	8.50	9.68	8.41	85.1	
P4104-23	WC-22A-EPH-3	34	1.12	8.66	9.78	9.21	93.4	
P4104-24	WC-22A-EPH-4	35	1.16	8.67	9.83	9.2	92.7	
P4104-25	WC-22A-EPH-5	36	1.16	8.70	9.86	9.24	92.9	
P4104-26	WC-23-EPH-1	37	1.17	8.56	9.73	8.72	88.2	
P4104-27	WC-23-EPH-2	38	1.15	8.83	9.98	9.00	88.9	
P4104-28	WC-23-EPH-3	39	1.15	8.81	9.96	8.65	85.1	
P4104-29	WC-23-EPH-4	40	1.15	8.37	9.52	8.56	88.5	
P4104-30	WC-23-EPH-5	41	1.16	8.61	9.77	8.63	86.8	
P4117-01	MOO-SOIL-PILE	42	1.18	8.68	9.86	8.86	88.5	
P4119-01	TRANSFOMER-OIL	43	1.00	1.00	2.00	2.00	100.0	oil sample
P4119-03	OILY-DEBRIS	44	1.00	1.00	2.00	2.00	100.0	oily-debris
P4120-01	T1	45	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4120-03	T2	46	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4120-05	SF-6-A	47	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4120-07	SF-6-B	48	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4120-09	COMP-1	49	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4120-11	COMP-2	50	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4120-13	RB24085	51	1.13	8.82	9.95	4.23	35.1	GEL MARTIX
P4121-01	COMP-1	52	1.13	8.51	9.64	7.12	70.4	
P4121-02	COMP-2	53	1.00	1.00	2.00	2.00	100.0	debris
P4121-03	ELZ-24-00090	54	1.00	1.00	2.00	2.00	100.0	oil sample
P4122-01	FDH013Y-1-1	55	1.00	1.00	2.00	2.00	100.0	oilc
P4122-02	FDH013Y-1-2	56	1.00	1.00	2.00	2.00	100.0	oilc



PERCENT SOLID

Supervisor: sohil
Analyst: jignesh
Date: 9/20/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 09/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 09/20/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132510

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4122-03	FDH013Y-2-1	57	1.00	1.00	2.00	2.00	100.0	oilc
P4122-04	FDH013Y-2-2	58	1.00	1.00	2.00	2.00	100.0	oilc
P4123-01	910-B	59	1.13	8.80	9.93	6.38	59.7	
P4123-02	613-D	60	1.00	1.00	2.00	2.00	100.0	debris
P4124-01	2049	61	1.00	1.00	2.00	2.00	100.0	wipe sample
P4124-02	2050	62	1.00	1.00	2.00	2.00	100.0	wipe sample
P4125-01	PL-01-09192024	63	1.15	8.68	9.83	9.53	96.5	
P4125-02	PL-01-09192024-E2	64	1.16	8.44	9.6	9.36	97.2	
P4126-01	NB-08-09192024	65	1.18	8.52	9.7	9.37	96.1	
P4126-02	NB-08-09192024-E2	66	1.11	8.85	9.96	9.27	92.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$