

Analysis Method: 9034

ANALYST: rubina

Parameter: Reactive Sulfide

SUPERVISOR REVIEW BY: Iwona

Run Number: LB132613

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3114
Starch Solution, 4L	W2977

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB163654BL		1	5.00	50	2.00	0.00	1.94	1.94	0.06	0.00	0.00	09/26/2024	10:20
2	P4144-02		1	5.06	50	2.00	0.00	1.90	1.90	0.10	0.04	3.16	09/26/2024	10:23
3	P4144-02DUP		1	5.06	50	2.00	0.00	1.90	1.90	0.10	0.04	3.16	09/26/2024	10:25
4	P4144-04		1	5.01	50	2.00	0.00	1.86	1.86	0.14	0.08	6.39	09/26/2024	10:28
5	P4148-02		1	5.04	50	2.00	0.00	1.86	1.86	0.14	0.08	6.35	09/26/2024	10:30
6	P4169-01		1	5.03	50	2.00	0.00	1.90	1.90	0.10	0.04	3.18	09/26/2024	10:33
7	P4172-03		1	5.02	50	2.00	0.00	1.88	1.88	0.12	0.06	4.78	09/26/2024	10:35
8	P4178-01		1	5.01	50	2.00	0.00	1.90	1.90	0.10	0.04	3.19	09/26/2024	10:38
9	P4178-02		1	5.04	50	2.00	0.00	1.86	1.86	0.14	0.08	6.35	09/26/2024	10:40
10	P4178-03		1	5.05	50	2.00	0.00	1.86	1.86	0.14	0.08	6.34	09/26/2024	10:42
11	P4190-04		1	5.03	50	2.00	0.00	1.88	1.88	0.12	0.06	4.77	09/26/2024	10:45
12	P4192-02		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.04	3.17	09/26/2024	10:48

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume

Analysis Method: 9034

Parameter: Reactive Sulfide

Run Number: LB132613

ANALYST: rubina

SUPERVISOR REVIEW BY: Iwona

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3114
Starch Solution, 4L	W2977

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
13	P4192-04		1	5.02	50	2.00	0.00	1.90	1.90	0.10	0.04	3.19	09/26/2024	10:50

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume