



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 10/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132772

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4300-05	SVOC-GPC-BLANK	1	1.00	1.00	2.00	2.00	100.0	
P4300-06	PEST-GPC-BLANK	2	1.00	1.00	2.00	2.00	100.0	
P4300-07	PEST-GPC-BLANK-SPIKE	3	1.00	1.00	2.00	2.00	100.0	
P4300-08	PCB-GPC-BLANK	4	1.00	1.00	2.00	2.00	100.0	
P4300-09	PCB-GPC-BLANK-SPIKE	5	1.00	1.00	2.00	2.00	100.0	
P4300-10	SVOC-GPC2-BLANK	6	1.00	1.00	2.00	2.00	100.0	
P4300-11	PEST-GPC2-BLANK	7	1.00	1.00	2.00	2.00	100.0	
P4300-12	PEST-GPC2-BLANK-SPIKE	8	1.00	1.00	2.00	2.00	100.0	
P4300-13	PCB-GPC2-BLANK	9	1.00	1.00	2.00	2.00	100.0	
P4300-14	PCB-GPC2-BLANK-SPIKE	10	1.00	1.00	2.00	2.00	100.0	
P4307-01	WS-05	11	1.00	1.00	2.00	2.00	100.0	wipe sample
P4307-02	WS-06	12	1.00	1.00	2.00	2.00	100.0	wipe sample
P4307-03	WS-07	13	1.00	1.00	2.00	2.00	100.0	wipe sample
P4307-04	WS-08	14	1.00	1.00	2.00	2.00	100.0	wipe sample
P4307-05	WS-09	15	1.00	1.00	2.00	2.00	100.0	wipe sample
P4307-06	WS-10	16	1.00	1.00	2.00	2.00	100.0	wipe sample
P4321-01	IB-6	17	1.17	8.40	9.57	6.1	58.7	
P4321-02	IB-6-EPH	18	1.18	8.71	9.89	5.92	54.4	
P4321-03	IB-6-VOC	19	1.18	8.51	9.69	8.62	87.4	
P4323-01	MOO-23-00323	20	1.00	1.00	2.00	2.00	100.0	debris
P4323-02	MOO-23-00324	21	1.15	8.81	9.96	8.71	85.8	
P4324-01	LAW-24-00169-00170	22	1.13	8.53	9.66	8.26	83.6	
P4324-02	LAW-24-00166	23	1.00	1.00	2.00	2.00	100.0	debris
P4324-03	JEC021I-1-1	24	1.00	1.00	2.00	2.00	100.0	oilc
P4324-04	JEC021I-1-2	25	1.00	1.00	2.00	2.00	100.0	oilc
P4326-01	VNJ-214	26	1.15	8.83	9.98	8.01	77.7	
P4327-01	JEH697J-1-1	27	1.00	1.00	2.00	2.00	100.0	oilc
P4327-02	JEH697J-1-2	28	1.00	1.00	2.00	2.00	100.0	oilc



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/04/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 10/05/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132772

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4327-03	KMH544K-1-1	29	1.00	1.00	2.00	2.00	100.0	oilc
P4327-04	KMH544K-1-2	30	1.00	1.00	2.00	2.00	100.0	oilc
P4328-01	4	31	1.00	1.00	2.00	2.00	100.0	wipe sample
P4328-02	5	32	1.00	1.00	2.00	2.00	100.0	wipe sample
P4328-03	6	33	1.00	1.00	2.00	2.00	100.0	wipe sample
P4328-04	7	34	1.00	1.00	2.00	2.00	100.0	wipe sample
P4328-05	8	35	1.00	1.00	2.00	2.00	100.0	wipe sample
P4328-06	9	36	1.00	1.00	2.00	2.00	100.0	wipe sample
P4328-07	CONCRETE-RO	37	1.00	1.00	2.00	2.00	100.0	CONCRETE sample

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$