

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/7/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 13:50
In Date: 10/06/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:50
Out Date: 10/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132786

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4294-01	MYDA19	1	1.18	8.44	9.62	9.57	99.4	
P4294-02	MYDA20	2	1.17	8.57	9.74	9.7	99.5	
P4294-03	MYDA21	3	1.18	8.74	9.92	9.82	98.9	
P4294-04	MYDA22	4	1.18	8.57	9.75	9.71	99.5	
P4294-05	MYDA23	5	1.17	8.60	9.77	9.74	99.7	
P4294-06	MYDA24	6	1.14	8.83	9.97	9.87	98.9	
P4294-07	MYDA25	7	1.15	8.80	9.95	9.9	99.4	
P4294-08	MYDA26	8	1.18	8.40	9.58	9.35	97.3	
P4294-09	MYDA27	9	1.18	8.51	9.69	9.58	98.7	
P4294-10	MYDA28	10	1.19	8.47	9.66	9.55	98.7	
P4294-11	MYDA29	11	1.15	8.83	9.98	9.9	99.1	
P4294-12	MYDA30	12	1.17	8.61	9.78	9.73	99.4	
P4294-13	MYDA31	13	1.15	8.65	9.8	9.72	99.1	
P4294-14	MYDA32	14	1.17	8.60	9.77	9.68	99.0	
P4294-15	MYDA33	15	1.18	8.44	9.62	9.55	99.2	
P4294-16	MYDA34	16	1.18	8.64	9.82	9.77	99.4	
P4294-17	MYDA35	17	1.16	8.77	9.93	9.8	98.5	
P4294-18	MYDA36	18	1.19	8.63	9.82	9.8	99.8	
P4294-19	MYDA37	19	1.18	8.53	9.71	9.66	99.4	
P4294-20	MYDA42	20	1.16	8.60	9.76	9.66	98.8	
P4294-21	MYDA42D	21	1.16	8.60	9.76	9.66	98.8	
P4294-22	MYDA42S	22	1.16	8.60	9.76	9.66	98.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$