



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/8/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 10/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:50
Out Date: 10/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132803

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4304-01	MYDAH3	1	1.16	8.46	9.62	9.37	97.0	
P4304-02	MYDAH4	2	1.16	8.51	9.67	9.54	98.5	
P4304-03	MYDAH5	3	1.12	8.42	9.54	9.38	98.1	
P4304-04	MYDAH6	4	1.15	8.75	9.9	9.7	97.7	
P4304-05	MYDAH7	5	1.15	8.53	9.68	9.64	99.5	
P4304-06	MYDAH8	6	1.16	8.46	9.62	9.58	99.5	
P4304-07	MYDAH9	7	1.13	8.79	9.92	9.87	99.4	
P4304-08	MYDAP2	8	1.15	8.64	9.79	9.66	98.5	
P4304-09	MYDAP2D	9	1.15	8.64	9.79	9.66	98.5	
P4304-10	MYDAP2S	10	1.15	8.64	9.79	9.66	98.5	
P4304-11	MYDAP3	11	1.13	8.43	9.56	9.44	98.6	
P4304-12	MYDAP4	12	1.14	8.84	9.98	9.92	99.3	
P4304-13	MYDAP5	13	1.14	8.75	9.89	9.85	99.5	
P4304-14	MYDAJ0	14	1.15	8.49	9.64	9.62	99.8	
P4304-15	MYDAJ1	15	1.12	8.61	9.73	9.71	99.8	
P4304-16	MYDAJ2	16	1.12	8.71	9.83	9.72	98.7	
P4304-17	MYDAJ3	17	1.14	8.73	9.87	9.84	99.7	
P4304-18	MYDAJ4	18	1.15	8.67	9.82	9.78	99.5	
P4304-19	MYDAJ6	19	1.15	8.52	9.67	9.62	99.4	
P4304-20	MYDAJ7	20	1.15	8.58	9.73	9.65	99.1	
P4304-21	MYDAJ8	21	1.13	8.47	9.6	9.55	99.4	
P4304-22	MYDAJ9	22	1.12	8.46	9.58	9.47	98.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$