



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/8/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 15:55
In Date: 10/07/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132809

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4308-01	MYDAN1	1	1.13	8.59	9.72	9.53	97.8	
P4308-02	MYDAN1D	2	1.13	8.59	9.72	9.53	97.8	
P4308-03	MYDAN1S	3	1.13	8.59	9.72	9.53	97.8	
P4308-04	MYDAP0	4	1.13	8.54	9.67	9.59	99.1	
P4308-05	MYDAP1	5	1.16	8.64	9.8	9.65	98.3	
P4308-06	MYDAP6	6	1.19	8.55	9.74	9.7	99.5	
P4308-07	MYDAP7	7	1.15	8.43	9.58	9.51	99.2	
P4308-08	MYDAP8	8	1.15	8.44	9.59	9.51	99.1	
P4308-09	MYDAP9	9	1.14	8.42	9.56	9.5	99.3	
P4308-10	MYDAQ0	10	1.13	8.68	9.81	9.75	99.3	
P4308-11	MYDAQ1	11	1.15	8.59	9.74	9.62	98.6	
P4308-12	MYDAQ2	12	1.17	8.35	9.52	9.32	97.6	
P4308-13	MYDAQ3	13	1.14	8.53	9.67	9.56	98.7	
P4308-14	MYDAQ4	14	1.14	8.83	9.97	9.82	98.3	
P4308-15	MYDAQ5	15	1.15	8.73	9.88	9.72	98.2	
P4308-16	MYDAQ7	16	1.16	8.57	9.73	9.62	98.7	
P4308-17	MYDAQ8	17	1.17	8.74	9.91	9.78	98.5	
P4308-18	MYDAQ9	18	1.15	8.58	9.73	9.63	98.8	
P4308-19	MYDAR0	19	1.15	8.38	9.53	9.33	97.6	
P4308-20	MYDAR1	20	1.16	8.58	9.74	9.65	99.0	
P4308-21	MYDAR2	21	1.17	8.44	9.61	9.44	98.0	
P4308-22	MYDAR3	22	1.17	8.75	9.92	9.63	96.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$