



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/10/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 12:25
In Date: 10/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 10/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132839

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4310-01	MYDAT3	1	1.14	8.58	9.72	9.45	96.9	
P4310-02	MYDAT4	2	1.15	8.37	9.52	9.38	98.3	
P4310-03	MYDAT5	3	1.18	8.50	9.68	9.51	98.0	
P4310-04	MYDAT5D	4	1.18	8.50	9.68	9.51	98.0	
P4310-05	MYDAT5S	5	1.18	8.50	9.68	9.51	98.0	
P4310-06	MYDAT6	6	1.17	8.55	9.72	9.55	98.0	
P4310-07	MYDAT7	7	1.13	8.80	9.93	9.68	97.2	
P4310-08	MYDAT8	8	1.18	8.60	9.78	9.61	98.0	
P4310-09	MYDAT9	9	1.18	8.62	9.8	9.7	98.8	
P4310-10	MYDAW0	10	1.18	8.47	9.65	9.53	98.6	
P4310-11	MYDAW1	11	1.18	8.70	9.88	9.72	98.2	
P4310-12	MYDAW2	12	1.17	8.59	9.76	9.57	97.8	
P4310-13	MYDAW3	13	1.16	8.61	9.77	9.57	97.7	
P4310-14	MYDAW4	14	1.18	8.55	9.73	9.6	98.5	
P4310-15	MYDAW5	15	1.16	8.40	9.56	9.37	97.7	
P4310-16	MYDAW6	16	1.17	8.49	9.66	9.53	98.5	
P4310-17	MYDAW7	17	1.18	8.63	9.81	9.71	98.8	
P4310-18	MYDAW8	18	1.14	8.40	9.54	9.34	97.6	
P4310-19	MYDAW9	19	1.18	8.45	9.63	9.5	98.5	
P4310-20	MYDAX0	20	1.15	8.73	9.88	9.77	98.7	
P4310-21	MYDAX1	21	1.18	8.44	9.62	9.3	96.2	
P4310-22	MYDAX2	22	1.18	8.75	9.93	9.78	98.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$