



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/10/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 12:45
In Date: 10/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 10/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132840

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4311-01	MYDAX3	1	1.13	8.66	9.79	9.56	97.3	
P4311-02	MYDAX4	2	1.18	8.49	9.67	9.55	98.6	
P4311-03	MYDAX5	3	1.16	8.48	9.64	9.51	98.5	
P4311-04	MYDAX6	4	1.15	8.63	9.78	9.53	97.1	
P4311-05	MYDAX7	5	1.15	8.43	9.58	9.32	96.9	
P4311-06	MYDAX8	6	1.15	8.35	9.5	9.27	97.2	
P4311-07	MYDAX9	7	1.15	8.83	9.98	9.72	97.1	
P4311-08	MYDAY0	8	1.19	8.61	9.8	9.5	96.5	
P4311-09	MYDAY1	9	1.15	8.65	9.8	9.52	96.8	
P4311-10	MYDAY2	10	1.15	8.65	9.8	9.35	94.8	
P4311-11	MYDAY3	11	1.12	8.70	9.82	9.5	96.3	
P4311-12	MYDAY4	12	1.15	8.44	9.59	9.34	97.0	
P4311-13	MYDAY5	13	1.16	8.61	9.77	9.63	98.4	
P4311-14	MYDAY6	14	1.15	8.36	9.51	9.21	96.4	
P4311-15	MYDAY7	15	1.17	8.60	9.77	9.69	99.1	
P4311-16	MYDAY8	16	1.17	8.56	9.73	9.57	98.1	
P4311-17	MYDAY9	17	1.15	8.67	9.82	9.75	99.2	
P4311-18	MYDAZ0	18	1.18	8.40	9.58	9.32	96.9	
P4311-19	MYDAZ1	19	1.18	8.78	9.96	9.68	96.8	
P4311-20	MYDB06	20	1.18	8.42	9.6	9.44	98.1	
P4311-21	MYDB06D	21	1.18	8.42	9.6	9.44	98.1	
P4311-22	MYDB06S	22	1.18	8.42	9.6	9.44	98.1	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$