

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 10/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:18
Out Date: 10/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132870

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4375-01	CONCRETE-PADS-1	1	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-03	CONCRETE-PADS-2	2	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-05	CONCRETE-PADS-3	3	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-07	CONCRETE-PADS-4	4	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-09	CONCRETE-PADS-5	5	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-11	CONCRETE-PADS-6	6	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-13	CONCRETE-PADS-7	7	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-15	CONCRETE-PADS-8	8	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-17	CONCRETE-PADS-9	9	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-19	CONCRETE-PADS-10	10	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-21	CONCRETE-PADS-11	11	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-23	FENCE-POST-1	12	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-25	FENCE-POST-2	13	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4375-27	FENCE-POST-3	14	1.00	1.00	2.00	2.00	100.0	CONCRETE sample
P4376-01	OR-02-100924	15	1.18	8.51	9.69	9.25	94.8	
P4376-02	OR-02-100924-E2	16	1.18	8.51	9.69	9.07	92.7	
P4377-01	HD-01-100924	17	1.17	8.53	9.7	8.9	90.6	
P4377-02	HD-01-100924-E2	18	1.15	8.83	9.98	9.03	89.2	
P4381-01	AT030P-SD04-100724-00	19	1.15	8.79	9.94	4.03	32.8	
P4381-02	AT030P-SD01-100724-00	20	1.15	8.70	9.85	8.44	83.8	
P4381-03	AT030P-SD02-100724-00	21	1.16	8.51	9.67	3.62	28.9	
P4383-01	34828	22	1.00	1.00	2.00	2.00	100.0	wipe sample
P4383-02	34829	23	1.00	1.00	2.00	2.00	100.0	wipe sample
P4383-03	34830	24	1.00	1.00	2.00	2.00	100.0	wipe sample
P4385-01	SP-1	25	1.17	8.60	9.77	9.15	92.8	
P4385-02	SP-1	26	1.14	8.48	9.62	9.06	93.4	
P4385-03	SP-2	27	1.15	8.41	9.56	9.09	94.4	
P4385-04	SP-2	28	1.18	8.39	9.57	9.15	95.0	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 10/10/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:18
Out Date: 10/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132870

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4385-05	SP-3	29	1.18	8.50	9.68	9.35	96.1	
P4385-06	SP-3	30	1.15	8.64	9.79	9.28	94.1	
P4385-07	SP-4	31	1.14	8.43	9.57	9.44	98.5	
P4385-08	SP-4	32	1.19	8.44	9.63	9.55	99.1	
P4385-09	SP-5	33	1.15	8.68	9.83	9.28	93.7	
P4385-10	SP-5	34	1.15	8.51	9.66	9.1	93.4	
P4385-11	SP-6	35	1.14	8.44	9.58	9.04	93.6	
P4385-12	SP-6	36	1.17	8.53	9.7	9.1	93.0	
P4385-13	SP-7	37	1.12	8.73	9.85	9.31	93.8	
P4385-14	SP-7	38	1.18	8.77	9.95	9.41	93.8	
P4385-15	SP-8	39	1.18	8.62	9.8	9.32	94.4	
P4385-16	SP-8	40	1.17	8.46	9.63	9.15	94.3	
P4385-17	SP-9	41	1.15	8.61	9.76	9.19	93.4	
P4385-18	SP-9	42	1.18	8.53	9.71	9.15	93.4	
P4385-19	SP-10	43	1.19	8.53	9.72	9.18	93.7	
P4385-20	SP-10	44	1.18	8.55	9.73	9.25	94.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$