



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/14/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 13:15
In Date: 10/11/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 10/12/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132888

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4350-01	EOAK8	1	1.15	8.84	9.99	4.76	40.8	
P4350-02	VHBLK001	2	1.00	1.00	2.00	2.00	100.0	vhblk
P4350-03	EOAK4	3	1.15	8.57	9.72	5.00	44.9	
P4350-04	EOAK5	4	1.15	8.40	9.55	7.87	80.0	
P4350-05	EOAK7	5	1.18	8.64	9.82	5.85	54.1	
P4350-06	EOAK3	6	1.16	8.72	9.88	7.17	68.9	
P4350-07	EOAK6	7	1.15	8.82	9.97	6.47	60.3	
P4350-08	EOAK9	8	1.14	8.47	9.61	6.71	65.8	
P4350-09	EOAJ9	9	1.13	8.74	9.87	8.18	80.7	
P4350-10	EOAK0	10	1.14	8.80	9.94	8.02	78.2	
P4350-11	EOAK2	11	1.18	8.63	9.81	6.67	63.6	
P4350-12	EOAK1	12	1.14	8.80	9.94	8.08	78.9	
P4350-13	EOAK1MS	13	1.14	8.80	9.94	8.08	78.9	
P4350-14	EOAK1MSD	14	1.14	8.80	9.94	8.08	78.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$