



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/17/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132959

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4412-01	KR-SS01	1	1.15	8.67	9.82	8.64	86.4	
P4412-02	KR-SS02	2	1.15	8.80	9.95	8.75	86.4	
P4412-03	KR-SS03	3	1.15	8.80	9.95	8.83	87.3	
P4412-04	KR-SS04	4	1.15	8.59	9.74	8.41	84.5	
P4412-05	KR-SS05	5	1.15	8.47	9.62	8.16	82.8	
P4412-06	KR-SS06	6	1.15	8.81	9.96	8.46	83.0	
P4412-07	KR-SS07	7	1.17	8.55	9.72	8.36	84.1	
P4412-08	KR-SS08	8	1.17	8.55	9.72	8.31	83.5	
P4412-09	KR-SS09	9	1.13	8.66	9.79	8.61	86.4	
P4412-10	KR-SS10	10	1.15	8.83	9.98	8.87	87.4	
P4412-11	KR-SS11	11	1.19	8.63	9.82	8.7	87.0	
P4412-12	KR-SS12	12	1.15	8.81	9.96	8.64	85.0	
P4412-13	KR-SS13	13	1.16	8.63	9.79	8.5	85.1	
P4412-14	KR-SS14	14	1.15	8.38	9.53	8.32	85.6	
P4412-15	KR-SS15	15	1.15	8.81	9.96	9.00	89.1	
P4412-16	KR-SS16	16	1.18	8.59	9.77	8.61	86.5	
P4412-17	KR-SS17	17	1.12	8.70	9.82	8.89	89.3	
P4412-18	KR-SS18	18	1.16	8.48	9.64	8.69	88.8	
P4413-01	KR-SS19	19	1.17	8.60	9.77	8.73	87.9	
P4413-02	KR-SS20	20	1.17	8.68	9.85	8.69	86.6	
P4413-03	KR-SS21	21	1.16	8.46	9.62	8.66	88.7	
P4413-04	KR-SS22	22	1.14	8.60	9.74	8.57	86.4	
P4413-05	KR-SS23	23	1.18	8.54	9.72	8.33	83.7	
P4413-06	KR-SS24	24	1.18	8.50	9.68	8.69	88.4	
P4413-07	KR-SS25	25	1.17	8.60	9.77	8.69	87.4	
P4413-08	KR-SS26	26	1.16	8.49	9.65	8.65	88.2	
P4413-09	KR-SS27	27	1.18	8.50	9.68	8.74	88.9	
P4413-10	KR-SS28	28	1.14	8.60	9.74	8.6	86.7	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/17/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 10/17/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB132959

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4413-11	KR-SS29	29	1.15	8.40	9.55	8.63	89.0	
P4413-12	KR-SS30	30	1.16	8.78	9.94	8.7	85.9	
P4413-13	KR-SS31	31	1.16	8.82	9.98	7.85	75.9	
P4413-14	KR-SS32	32	1.16	8.63	9.79	8.68	87.1	
P4413-15	KR-SS33	33	1.13	8.81	9.94	8.79	86.9	
P4417-01	HOPPER-A-SOLID	34	1.14	8.61	9.75	8.38	84.1	
P4417-04	HOPPER-B-SOLID	35	1.11	8.86	9.97	9.22	91.5	
P4417-05	HOPPER-B-SOLID	36	1.15	8.81	9.96	9.05	89.7	
P4417-07	HOPPER-C-SOLID	37	1.15	8.78	9.93	8.06	78.7	
P4417-08	HOPPER-C-SOLID	38	1.15	8.80	9.95	8.48	83.3	
P4419-01	IB-1	39	1.16	8.40	9.56	7.04	70.0	
P4419-02	IB-1-EPH	40	1.14	8.61	9.75	7.2	70.4	
P4419-03	IB-1-VOC	41	1.11	8.53	9.64	6.98	68.8	
P4421-01	EO-02-101624	42	1.15	8.58	9.73	8.88	90.1	
P4421-02	EO-02-101624-E2	43	1.15	8.83	9.98	9.09	89.9	
P4422-01	EO-01-101624	44	1.16	8.50	9.66	9.22	94.8	
P4422-02	EO-01-101624-E2	45	1.17	8.73	9.9	9.42	94.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$