



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/23/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 10/22/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:49
Out Date: 10/23/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133060

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4477-01	MYE450	1	1.15	8.53	9.68	8.98	91.8	
P4477-02	MYE451	2	1.15	8.70	9.85	9.65	97.7	
P4477-03	MYE452	3	1.16	8.52	9.68	9.05	92.6	
P4477-04	MYE453	4	1.17	8.54	9.71	9.47	97.2	
P4477-05	MYE453D	5	1.17	8.54	9.71	9.47	97.2	
P4477-06	MYE453S	6	1.17	8.54	9.71	9.47	97.2	
P4477-07	MYE454	7	1.15	8.80	9.95	9.27	92.3	
P4477-08	MYE455	8	1.14	8.48	9.62	8.94	92.0	
P4477-09	MYE456	9	1.15	8.68	9.83	9.59	97.2	
P4477-10	MYE457	10	1.18	8.63	9.81	9.6	97.6	
P4477-11	MYE458	11	1.14	8.53	9.67	9.44	97.3	
P4477-12	MYE459	12	1.15	8.44	9.59	9.32	96.8	
P4477-13	MYE460	13	1.15	8.62	9.77	9.52	97.1	
P4477-14	MYE461	14	1.19	8.50	9.69	9.6	98.9	
P4477-15	MYE462	15	1.18	8.49	9.67	9.2	94.5	
P4477-16	MYE463	16	1.15	8.43	9.58	8.85	91.3	
P4477-17	MYE464	17	1.16	8.69	9.85	9.14	91.8	
P4477-18	MYE465	18	1.17	8.78	9.95	9.5	94.9	
P4477-19	MYE466	19	1.17	8.71	9.88	9.63	97.1	
P4477-20	MYE467	20	1.16	8.58	9.74	9.52	97.4	
P4477-21	MYE468	21	1.17	8.70	9.87	9.4	94.6	
P4477-22	MYE469	22	1.18	8.75	9.93	9.47	94.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$