



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/24/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 12:30
In Date: 10/23/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:37
Out Date: 10/24/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133066

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4479-01	MYE470	1	1.15	8.53	9.68	9.34	96.0	
P4479-02	MYE471	2	1.15	8.52	9.67	9.29	95.5	
P4479-03	MYE472	3	1.15	8.50	9.65	9.18	94.5	
P4479-04	MYE473	4	1.17	8.53	9.7	9.37	96.1	
P4479-05	MYE474	5	1.17	8.62	9.79	9.39	95.4	
P4479-06	MYE475	6	1.16	8.47	9.63	9.08	93.5	
P4479-07	MYE476	7	1.17	8.56	9.73	9.33	95.3	
P4479-08	MYE477	8	1.17	8.48	9.65	9.33	96.2	
P4479-09	MYE478	9	1.16	8.40	9.56	9.24	96.2	
P4479-10	MYE479	10	1.17	8.53	9.7	9.33	95.7	
P4479-11	MYE480	11	1.18	8.57	9.75	9.44	96.4	
P4479-12	MYE481	12	1.17	8.49	9.66	9.26	95.3	
P4479-13	MYE481D	13	1.17	8.49	9.66	9.26	95.3	
P4479-14	MYE481S	14	1.17	8.49	9.66	9.26	95.3	
P4479-15	MYE482	15	1.15	8.65	9.8	9.65	98.3	
P4479-16	MYE483	16	1.15	8.66	9.81	9.57	97.2	
P4479-17	MYE484	17	1.14	8.65	9.79	9.43	95.8	
P4479-18	MYE485	18	1.16	8.61	9.77	9.55	97.4	
P4479-19	MYE486	19	1.16	8.60	9.76	9.16	93.0	
P4479-20	MYE487	20	1.17	8.73	9.9	9.53	95.8	
P4479-21	MYE488	21	1.16	8.62	9.78	9.36	95.1	
P4479-22	MYE489	22	1.17	8.73	9.9	9.5	95.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$