



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/24/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 13:20
In Date: 10/23/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:50
Out Date: 10/24/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133069

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4435-01	BHK68	1	1.15	8.84	9.99	9.26	91.7	
P4435-02	BHK68MS	2	1.15	8.84	9.99	9.26	91.7	
P4435-03	BHK68MSD	3	1.15	8.84	9.99	9.26	91.7	
P4481-01	MYE4B0	4	1.15	8.72	9.87	9.7	98.1	
P4481-02	MYE4B1	5	1.18	8.43	9.61	9.36	97.0	
P4481-03	MYE4B2	6	1.15	8.67	9.82	9.47	96.0	
P4481-04	MYE4B3	7	1.17	8.58	9.75	9.48	96.9	
P4481-05	MYE4B4	8	1.16	8.46	9.62	9.35	96.8	
P4481-06	MYE4B5	9	1.15	8.72	9.87	9.63	97.2	
P4481-07	MYE4B6	10	1.14	8.63	9.77	9.35	95.1	
P4481-08	MYE4B7	11	1.17	8.69	9.86	9.57	96.7	
P4481-09	MYE4B8	12	1.15	8.56	9.71	9.55	98.1	
P4481-10	MYE4B9	13	1.17	8.38	9.55	9.3	97.0	
P4481-11	MYE4C0	14	1.15	8.71	9.86	9.57	96.7	
P4481-12	MYE4C1	15	1.17	8.53	9.7	9.27	95.0	
P4481-13	MYE4C1D	16	1.17	8.53	9.7	9.27	95.0	
P4481-14	MYE4C1S	17	1.17	8.53	9.7	9.27	95.0	
P4481-15	MYE4C2	18	1.16	8.36	9.52	9.16	95.7	
P4481-16	MYE4C3	19	1.17	8.62	9.79	9.56	97.3	
P4481-17	MYE4C4	20	1.17	8.41	9.58	9.15	94.9	
P4481-18	MYE4C5	21	1.16	8.56	9.72	9.4	96.3	
P4481-19	MYE4C6	22	1.16	8.37	9.53	9.34	97.7	
P4481-20	MYE4C7	23	1.17	8.61	9.78	9.4	95.6	
P4481-21	MYE4C8	24	1.16	8.71	9.87	9.57	96.6	
P4481-22	MYE4C9	25	1.17	8.63	9.8	9.59	97.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$