

Analysis Method: 9034

ANALYST: rubina

Parameter: Reactive Sulfide

SUPERVISOR REVIEW BY: Iwona

Run Number: LB133098

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3114
Starch Solution, 4L	W3149

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB164329BL		1	5.00	50	2.00	0.00	1.92	1.92	0.08	0.00	0.00	10/24/2024	11:20
2	P4485-02		1	5.06	50	2.00	0.00	1.90	1.90	0.10	0.02	1.58	10/24/2024	11:23
3	P4485-02DUP		1	5.06	50	2.00	0.00	1.90	1.90	0.10	0.02	1.58	10/24/2024	11:25
4	P4487-04		1	5.02	50	2.00	0.00	1.86	1.86	0.14	0.06	4.78	10/24/2024	11:28
5	P4487-08		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	10/24/2024	11:30
6	P4508-04		1	5.03	50	2.00	0.00	1.88	1.88	0.12	0.04	3.18	10/24/2024	11:33
7	P4508-08		1	5.02	50	2.00	0.00	1.86	1.86	0.14	0.06	4.78	10/24/2024	11:36
8	P4508-12		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	10/24/2024	11:38
9	P4511-02		1	5.07	50	2.00	0.00	1.88	1.88	0.12	0.04	3.16	10/24/2024	11:40
10	P4513-01		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.02	1.58	10/24/2024	11:43
11	P4513-02		1	5.03	50	2.00	0.00	1.86	1.86	0.14	0.06	4.77	10/24/2024	11:46
12	P4513-03		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.02	1.58	10/24/2024	11:48

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume

Analysis Method: 9034

Parameter: Reactive Sulfide

Run Number: LB133098

ANALYST: rubina

SUPERVISOR REVIEW BY: Iwona

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3114
Starch Solution, 4L	W3149

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
13	P4515-02		1	5.02	50	2.00	0.00	1.84	1.84	0.16	0.08	6.37	10/24/2024	11:50
14	P4516-01		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	10/24/2024	11:53

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume