



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/25/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 10/24/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 10/25/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133108

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4509-01	AU-06-10232024	1	1.15	8.65	9.8	9.02	91.0	
P4518-01	EX-5-TPH-1	2	1.18	8.55	9.73	7.89	78.5	
P4518-02	EX-5-TPH-2	3	1.15	8.84	9.99	8.15	79.2	
P4518-03	EX-5-TPH-3	4	1.19	8.61	9.8	8.43	84.1	
P4518-04	EX-5-TPH-4	5	1.18	8.53	9.71	8.02	80.2	
P4518-05	EX-5-TPH-5	6	1.16	8.51	9.67	7.45	73.9	
P4531-01	OR-03-102424	7	1.19	8.61	9.8	9.31	94.3	
P4531-02	OR-03-102424-E2	8	1.19	8.63	9.82	9.34	94.4	
P4545-01	VNJ-215	15	1.13	8.69	9.82	8.69	87.0	
P4545-02	VNJ-215	16	1.13	8.69	9.82	8.69	87.0	
P4547-01	BP-F-21	9	1.14	8.82	9.96	8.98	88.9	
P4547-02	BP-F-21-EPH	10	1.18	8.66	9.84	8.83	88.3	
P4547-03	BP-F-21-VOC	11	1.13	8.75	9.88	8.81	87.8	
P4547-05	BP-F-20	12	1.16	8.80	9.96	8.78	86.6	
P4547-06	BP-F-20-EPH	13	1.18	8.43	9.61	8.45	86.2	
P4547-07	BP-F-20-VOC	14	1.14	8.85	9.99	8.82	86.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$