



PERCENT SOLID

Supervisor: Iwona  
Analyst: jignesh  
Date: 10/28/2024

OVENTEMP IN Celsius(°C): 107  
Time IN: 18:15  
In Date: 10/25/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:38  
Out Date: 10/26/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: % SOLID- OVEN

QC:LB133131

| Lab ID   | Client SampleID       | Dish # | Dish Wt (g) (A) | Sample Wt (g) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments    |
|----------|-----------------------|--------|-----------------|---------------|--------------------------|----------------------------|---------|-------------|
| P4561-01 | BP-F-19               | 1      | 1.16            | 8.48          | 9.64                     | 8.73                       | 89.3    |             |
| P4561-02 | BP-F-19-EPH           | 2      | 1.16            | 8.50          | 9.66                     | 8.77                       | 89.5    |             |
| P4561-03 | BP-F-19-VOC           | 3      | 1.15            | 8.80          | 9.95                     | 9.16                       | 91.0    |             |
| P4561-05 | BP-F-18               | 4      | 1.15            | 8.40          | 9.55                     | 8.36                       | 85.8    |             |
| P4561-06 | BP-F-18-EPH           | 5      | 1.13            | 8.44          | 9.57                     | 8.29                       | 84.8    |             |
| P4561-07 | BP-F-18-VOC           | 6      | 1.12            | 8.71          | 9.83                     | 8.52                       | 85.0    |             |
| P4562-05 | SVOC-GPC-BLANK        | 7      | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-06 | PEST-GPC-BLANK        | 8      | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-07 | PEST-GPC-BLANK-SPIKE  | 9      | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-08 | PCB-GPC-BLANK         | 10     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-09 | PCB-GPC-BLANK-SPIKE   | 11     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-10 | SVOC-GPC2-BLANK       | 12     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-11 | PEST-GPC2-BLANK       | 13     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-12 | PEST-GPC2-BLANK-SPIKE | 14     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-13 | PCB-GPC2-BLANK        | 15     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4562-14 | PCB-GPC2-BLANK-SPIKE  | 16     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   |             |
| P4563-01 | 100724-A              | 26     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-02 | 100724-B              | 27     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-03 | 100724-C              | 28     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-04 | 100724-D              | 29     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-05 | 100824-A              | 30     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-06 | 100924-A              | 31     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-07 | 100924-B              | 32     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-08 | 101424-A              | 33     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-09 | 102324-A              | 34     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-10 | 102524-A              | 35     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-11 | FITTING-1             | 36     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |
| P4563-12 | FITTING-2             | 37     | 1.00            | 1.00          | 2.00                     | 2.00                       | 100.0   | wipe sample |



PERCENT SOLID

Supervisor: Iwona  
Analyst: jignesh  
Date: 10/28/2024

OVENTEMP IN Celsius(°C): 107  
Time IN: 18:15  
In Date: 10/25/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:38  
Out Date: 10/26/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: % SOLID- OVEN

QC:LB133131

| Lab ID   | Client SampleID | Dish # | Dish Wt(g) (A) | Sample Wt(g) | Dish + Sample Wt(g) (B) | Dish+Dry Sample Wt(g) (C) | % Solid | Comments                   |
|----------|-----------------|--------|----------------|--------------|-------------------------|---------------------------|---------|----------------------------|
| P4565-01 | 2054            | 38     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | wipe sample                |
| P4565-02 | 2055            | 39     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | wipe sample                |
| P4565-03 | 2056            | 40     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | wipe sample                |
| P4565-04 | 2057            | 41     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | wipe sample                |
| P4565-05 | 2058            | 42     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | wipe sample                |
| P4565-06 | 2059            | 43     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | wipe sample                |
| P4566-01 | HD-01-102524    | 44     | 1.14           | 8.57         | 9.71                    | 9.47                      | 97.2    |                            |
| P4566-02 | HD-01-102524-E2 | 45     | 1.13           | 8.84         | 9.97                    | 9.42                      | 93.8    |                            |
| P4567-01 | WC-1            | 17     | 1.16           | 8.50         | 9.66                    | 8.89                      | 90.9    |                            |
| P4567-02 | WC-1-EPH        | 18     | 1.13           | 8.67         | 9.8                     | 8.99                      | 90.7    |                            |
| P4567-03 | WC-1-VOC        | 19     | 1.18           | 8.35         | 9.53                    | 9.05                      | 94.3    |                            |
| P4567-05 | WC-2            | 20     | 1.11           | 8.58         | 9.69                    | 8.99                      | 91.8    |                            |
| P4567-06 | WC-2-EPH        | 21     | 1.18           | 8.79         | 9.97                    | 9.12                      | 90.3    |                            |
| P4567-07 | WC-2-VOC        | 22     | 1.16           | 8.46         | 9.62                    | 9.46                      | 98.1    |                            |
| P4567-09 | WC-3            | 23     | 1.17           | 8.41         | 9.58                    | 8.69                      | 89.4    |                            |
| P4567-10 | WC-3-EPH        | 24     | 1.19           | 8.78         | 9.97                    | 9.17                      | 90.9    |                            |
| P4567-11 | WC-3-VOC        | 25     | 1.14           | 8.82         | 9.96                    | 9.19                      | 91.3    |                            |
| P4574-01 | GRAVEL-1        | 50     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | stone sample, 100 % solids |
| P4574-02 | GRAVEL-1-E2     | 51     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | stone sample, 100 % solids |
| P4574-04 | GRAVEL-2        | 52     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | stone sample, 100 % solids |
| P4574-05 | GRAVEL-2-E2     | 53     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | stone sample, 100 % solids |
| P4574-07 | CONCRETE-1      | 54     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | CONCRETE sample            |
| P4574-09 | CONCRETE-2      | 55     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | CONCRETE sample            |
| P4574-11 | CONCRETE-3      | 56     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | CONCRETE sample            |
| P4575-01 | PL-02-102424    | 57     | 1.15           | 8.42         | 9.57                    | 8.84                      | 91.3    |                            |
| P4575-02 | PL-02-102424-E2 | 58     | 1.15           | 8.83         | 9.98                    | 8.99                      | 88.8    |                            |
| P4576-01 | 1097            | 59     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | WIPE SAMPLE                |
| P4576-02 | 1098            | 60     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | WIPE Sample                |



PERCENT SOLID

Supervisor: Iwona  
Analyst: jignesh  
Date: 10/28/2024

OVENTEMP IN Celsius(°C): 107  
Time IN: 18:15  
In Date: 10/25/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:38  
Out Date: 10/26/2024  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: % SOLID- OVEN

QC:LB133131

| Lab ID   | Client SampleID | Dish # | Dish Wt (g) (A) | Sample Wt (g) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments |
|----------|-----------------|--------|-----------------|---------------|--------------------------|----------------------------|---------|----------|
| P4577-01 | TR-05-102524    | 61     | 1.18            | 8.40          | 9.58                     | 9.03                       | 93.5    |          |
| P4577-02 | TR-05-102524-E2 | 62     | 1.17            | 8.53          | 9.7                      | 9.03                       | 92.1    |          |
| P4578-03 | WB-305-TOP      | 46     | 1.15            | 8.81          | 9.96                     | 5.89                       | 53.8    |          |
| P4578-04 | WB-305-TOP-1    | 47     | 1.12            | 8.82          | 9.94                     | 6.3                        | 58.7    |          |
| P4578-05 | WB-305-BOT      | 48     | 1.15            | 8.81          | 9.96                     | 8.38                       | 82.1    |          |
| P4578-07 | WB-305-BOT-1    | 49     | 1.13            | 8.86          | 9.99                     | 9.49                       | 94.4    |          |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$