



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/28/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 13:05
In Date: 10/26/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 10/27/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133144

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4523-01	MYE4T2	1	1.16	8.78	9.94	9.72	97.5	
P4523-02	MYE4T3	2	1.15	8.41	9.56	9.5	99.3	
P4523-03	MYE4T4	3	1.15	8.63	9.78	9.34	94.9	
P4523-04	MYE4T5	4	1.16	8.42	9.58	9.39	97.7	
P4523-05	MYE4T6	5	1.16	8.47	9.63	9.45	97.9	
P4523-06	MYE4T7	6	1.16	8.70	9.86	9.72	98.4	
P4523-07	MYE4T8	7	1.18	8.37	9.55	9.46	98.9	
P4523-08	MYE4T8D	8	1.18	8.37	9.55	9.46	98.9	
P4523-09	MYE4T8S	9	1.18	8.37	9.55	9.46	98.9	
P4523-10	MYE4T9	10	1.18	8.66	9.84	9.66	97.9	
P4523-11	MYE4W8	11	1.16	8.67	9.83	9.67	98.2	
P4523-12	MYE4W9	12	1.16	8.65	9.81	9.63	97.9	
P4523-13	MYE4X0	13	1.16	8.48	9.64	9.53	98.7	
P4523-14	MYE4X1	14	1.18	8.73	9.91	9.66	97.1	
P4523-15	MYE4X2	15	1.15	8.80	9.95	9.61	96.1	
P4523-16	MYE4X3	16	1.18	8.64	9.82	9.6	97.5	
P4523-17	MYE4X4	17	1.17	8.59	9.76	9.53	97.3	
P4523-18	MYE4X5	18	1.16	8.59	9.75	9.61	98.4	
P4523-19	MYE4X6	19	1.17	8.70	9.87	9.55	96.3	
P4523-20	MYE4X7	20	1.18	8.76	9.94	9.77	98.1	
P4523-21	MYE4X8	21	1.15	8.57	9.72	9.53	97.8	
P4523-22	MYE4X9	22	1.15	8.43	9.58	9.32	96.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$