



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/28/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 14:55
In Date: 10/27/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:38
Out Date: 10/28/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133155

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4570-01	MYE531	1	1.14	8.58	9.72	9.38	96.0	
P4570-02	MYE531D	2	1.14	8.58	9.72	9.38	96.0	
P4570-03	MYE531S	3	1.14	8.58	9.72	9.38	96.0	
P4570-04	MYE537	4	1.18	8.49	9.67	9.31	95.8	
P4570-05	MYE538	5	1.15	8.46	9.61	8.85	91.0	
P4570-06	MYE539	6	1.15	8.54	9.69	7.61	75.6	
P4570-07	MYE540	7	1.15	8.77	9.92	9.73	97.8	
P4570-08	MYE541	8	1.15	8.53	9.68	9.62	99.3	
P4570-09	MYE542	9	1.15	8.63	9.78	9.45	96.2	
P4570-10	MYE563	10	1.18	8.44	9.62	9.4	97.4	
P4570-11	MYE564	11	1.18	8.69	9.87	9.63	97.2	
P4570-12	MYE565	12	1.18	8.66	9.84	9.65	97.8	
P4570-13	MYE566	13	1.12	8.77	9.89	9.73	98.2	
P4570-14	MYE567	14	1.12	8.59	9.71	9.48	97.3	
P4570-15	MYE568	15	1.18	8.74	9.92	9.63	96.7	
P4570-16	MYE569	16	1.18	8.40	9.58	9.37	97.5	
P4570-17	MYE570	17	1.15	8.69	9.84	9.55	96.7	
P4570-18	MYE571	18	1.15	8.80	9.95	9.68	96.9	
P4570-19	MYE572	19	1.18	8.53	9.71	9.51	97.7	
P4570-20	MYE573	20	1.18	8.75	9.93	9.66	96.9	
P4570-21	MYE574	21	1.12	8.70	9.82	9.35	94.6	
P4570-22	MYE575	22	1.13	8.61	9.74	9.36	95.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$