



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/30/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:25
In Date: 10/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/30/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB133200

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4609-01	COPI0	1	1.14	8.82	9.96	7.38	70.7	
P4609-02	COPI2	2	1.18	8.42	9.6	9.15	94.7	
P4609-03	COPI6	3	1.14	8.64	9.78	8.53	85.5	
P4609-04	COPI8	4	1.15	8.84	9.99	8.31	81.0	
P4609-05	COPI9	5	1.18	8.63	9.81	8.68	86.9	
P4609-06	CC0P1	6	1.13	8.70	9.83	9.42	95.3	
P4609-07	CC0P3	7	1.13	8.61	9.74	8.68	87.7	
P4609-08	CC0P5	8	1.17	8.55	9.72	8.77	88.9	
P4609-09	CC0P7	9	1.18	8.33	9.51	8.51	88.0	
P4609-10	CC0P9	10	1.15	8.56	9.71	8.71	88.3	
P4609-11	CC0Q1	11	1.17	8.53	9.7	8.68	88.0	
P4609-12	CC0Q6	12	1.16	8.69	9.85	9.22	92.8	
P4609-13	CC0Q8	13	1.16	8.50	9.66	8.35	84.6	
P4609-14	CC0R3	14	1.17	8.56	9.73	7.84	77.9	
P4609-15	CC0R4	15	1.16	8.42	9.58	7.19	71.6	
P4609-16	CC0R5	16	1.18	8.77	9.95	9.55	95.4	
P4609-17	CC0R7	17	1.15	8.37	9.52	7.95	81.2	
P4609-18	P4609-17MS	18	1.15	8.37	9.52	7.95	81.2	
P4609-19	P4609-17MSD	19	1.15	8.37	9.52	7.95	81.2	
P4610-01	COPI7	20	1.17	8.65	9.82	9.06	91.2	
P4610-02	CC0P0	21	1.15	8.57	9.72	9.1	92.8	
P4610-03	P4610-02MS	22	1.15	8.57	9.72	9.1	92.8	
P4610-04	P4610-02MSD	23	1.15	8.57	9.72	9.1	92.8	
P4610-05	CC0P2	24	1.15	8.78	9.93	9.58	96.0	
P4610-06	CC0P4	25	1.14	8.41	9.55	9.17	95.5	
P4610-07	CC0P6	26	1.14	8.64	9.78	9.11	92.2	
P4610-08	CC0P8	27	1.18	8.57	9.75	8.71	87.9	
P4610-09	CC0Q0	28	1.16	8.40	9.56	8.7	89.8	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 10/30/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:25
In Date: 10/29/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 10/30/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB133200

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4610-10	CC0Q2	29	1.13	8.54	9.67	8.43	85.5	
P4610-11	CC0Q3	30	1.13	8.47	9.6	6.8	66.9	
P4610-12	CC0Q4	31	1.16	8.42	9.58	7.84	79.3	
P4610-13	CC0Q5	32	1.15	8.38	9.53	7.64	77.4	
P4610-14	CC0Q7	33	1.15	8.78	9.93	8.61	85.0	
P4610-15	CC0Q9	34	1.14	8.79	9.93	9.57	95.9	
P4610-16	CC0R0	35	1.13	8.83	9.96	8.82	87.1	
P4610-17	CC0R1	36	1.15	8.66	9.81	7.9	77.9	
P4610-18	CC0R2	37	1.15	8.64	9.79	7.53	73.8	
P4610-19	CC0R8	38	1.15	8.41	9.56	8.85	91.6	
P4610-20	CC0R6	39	1.15	8.62	9.77	9.3	94.5	
P4610-21	P4610-20MS	40	1.15	8.62	9.77	9.3	94.5	
P4610-22	P4610-20MSD	41	1.15	8.62	9.77	9.3	94.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$