



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 11/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133361

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4716-01	A0BE3	1	1.15	8.83	9.98	8.14	79.2	
P4716-02	A0BE4	2	1.15	8.41	9.56	8.06	82.2	
P4716-03	A0BE5	3	1.17	8.62	9.79	6.03	56.4	
P4716-04	A0BE5MS	4	1.17	8.62	9.79	6.03	56.4	
P4716-05	A0BE5MSD	5	1.17	8.62	9.79	6.03	56.4	
P4771-05	P-1	6	1.12	8.74	9.86	9.37	94.4	
P4771-07	P-2	7	1.12	8.70	9.82	9.41	95.3	
P4779-01	SOIL-01	8	1.16	8.45	9.61	9.34	96.8	
P4780-05	SVOC-GPC-BLANK	9	1.00	1.00	2.00	2.00	100.0	
P4780-06	PEST-GPC-BLANK	10	1.00	1.00	2.00	2.00	100.0	
P4780-07	PEST-GPC-BLANK-SPIKE	11	1.00	1.00	2.00	2.00	100.0	
P4780-08	PCB-GPC-BLANK	12	1.00	1.00	2.00	2.00	100.0	
P4780-09	PCB-GPC-BLANK-SPIKE	13	1.00	1.00	2.00	2.00	100.0	
P4780-10	SVOC-GPC2-BLANK	14	1.00	1.00	2.00	2.00	100.0	
P4780-11	PEST-GPC2-BLANK	15	1.00	1.00	2.00	2.00	100.0	
P4780-12	PEST-GPC2-BLANK-SPIKE	16	1.00	1.00	2.00	2.00	100.0	
P4780-13	PCB-GPC2-BLANK	17	1.00	1.00	2.00	2.00	100.0	
P4780-14	PCB-GPC2-BLANK-SPIKE	18	1.00	1.00	2.00	2.00	100.0	
P4788-01	BP-G3	19	1.15	8.82	9.97	9.18	91.0	
P4788-02	BP-G3-EPH	20	1.13	8.68	9.81	8.85	88.9	
P4788-03	BP-G3-VOC	21	1.18	8.49	9.67	8.76	89.3	
P4789-01	BP-F21	22	1.13	8.76	9.89	8.98	89.6	
P4789-02	BP-F21-EPH	23	1.17	8.58	9.75	8.74	88.2	
P4789-03	BP-F21-VOC	24	1.17	8.58	9.75	8.74	88.2	
P4791-01	ARS20-0016	25	1.00	1.00	2.00	2.00	100.0	debris
P4792-02	DIESLE-DRUM	26	1.00	1.00	2.00	2.00	100.0	debris
P4792-03	OILY-STONE-DRUM	27	1.00	1.00	2.00	2.00	100.0	oily stone - drum
P4793-01	M00-24-00345	28	1.00	1.00	2.00	2.00	100.0	debris



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/11/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 11/08/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 11/09/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133361

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4794-01	KZA443F-END-1-1	29	1.00	1.00	2.00	2.00	100.0	pilc
P4794-02	KZA443F-END-1-2	30	1.00	1.00	2.00	2.00	100.0	pilc
P4795-01	LAW-23-00189	31	1.14	8.59	9.73	9.06	92.2	
P4796-01	TP-1	32	1.18	8.60	9.78	9.34	94.9	
P4796-02	TP-1	33	1.16	8.49	9.65	9.25	95.3	
P4796-03	TP-2	34	1.13	8.70	9.83	9.47	95.9	
P4796-04	TP-2	35	1.19	8.62	9.81	9.44	95.7	
P4796-05	TP-3	36	1.15	8.80	9.95	9.47	94.5	
P4796-06	TP-3	37	1.16	8.56	9.72	9.3	95.1	
P4796-07	TP-4	38	1.17	8.57	9.74	9.43	96.4	
P4796-08	TP-4	39	1.15	8.82	9.97	9.7	96.9	
P4798-01	MH-6	40	1.14	8.81	9.95	9.00	89.2	
P4798-02	MH-6-EPH	41	1.13	8.74	9.87	8.82	88.0	
P4798-03	MH-6-VOC	42	1.16	8.50	9.66	8.66	88.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$