



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/14/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 11/13/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 11/14/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133426

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4810-08	6	12	1.15	8.80	9.95	8.53	83.9	
P4824-01	T-20	1	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4824-02	T-20-E2	2	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4825-01	STOCK-PILE	3	1.15	8.46	9.61	8.73	89.6	
P4825-02	STOCK-PILE-E2	4	1.16	8.70	9.86	8.96	89.7	
P4826-01	GAS-BUR-2005	5	1.00	1.00	2.00	2.00	100.0	debris
P4833-01	MH-731	6	1.15	8.77	9.92	8.72	86.3	
P4833-02	MH-731-EPH	7	1.17	8.71	9.88	8.66	86.0	
P4833-03	MH-731-VOC	8	1.15	8.82	9.97	8.81	86.8	
P4833-05	TP-2	9	1.15	8.40	9.55	9.06	94.2	
P4833-06	TP-2-EPH	10	1.17	8.50	9.67	9.11	93.4	
P4833-07	TP-2-VOC	11	1.16	8.50	9.66	9.17	94.2	
P4834-03	SLP-1-STONE	13	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4834-04	SLP-1-STONE-E2	14	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4835-03	345-2-STONE	15	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4835-04	345-2-STONE-E2	16	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4836-03	T3-STONE	17	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4836-04	T3-STONE-E2	18	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4837-03	3TRK-STONE	19	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4837-04	3TRK-STONE-E2	20	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4838-03	T1-STONE	21	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids
P4838-04	T1-STONE-E2	22	1.00	1.00	2.00	2.00	100.0	stone sample, 100 % solids

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$