

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 11/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 11/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133466

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4851-01	OK-02-11142024	1	1.15	8.40	9.55	9.06	94.2	
P4851-02	OK-02-11142024-E2	2	1.16	8.50	9.66	9.32	96.0	
P4860-01	DUP-01	15	1.15	8.38	9.53	9.26	96.8	
P4860-02	PH2-BPT-001	16	1.18	8.50	9.68	9.41	96.8	
P4860-03	PH2-BPT-002	17	1.16	8.47	9.63	9.34	96.6	
P4860-04	PH2-BPT-003	18	1.15	8.49	9.64	9.45	97.8	
P4860-05	PH2-BPT-004	19	1.16	8.78	9.94	9.77	98.1	
P4860-06	PH2-BPT-009	20	1.19	8.59	9.78	9.63	98.3	
P4860-07	PH2-BPT-008	21	1.16	8.48	9.64	9.37	96.8	
P4860-08	PH2-BPT-007	22	1.16	8.48	9.64	9.47	98.0	
P4860-09	PH2-BPT-006	23	1.15	8.81	9.96	9.74	97.5	
P4860-10	PH2-BPT-005	24	1.13	8.85	9.98	9.77	97.6	
P4862-01	Y2404-0033-END-1-1	45	1.00	1.00	2.00	2.00	100.0	oilc
P4862-02	Y2404-0033-END-1-2	46	1.00	1.00	2.00	2.00	100.0	oilc
P4862-03	Y2404-END-1-1	47	1.00	1.00	2.00	2.00	100.0	oilc
P4862-04	Y2404-END-1-2	48	1.00	1.00	2.00	2.00	100.0	oilc
P4862-05	KMH4671-END-1-1	49	1.00	1.00	2.00	2.00	100.0	oilc
P4862-06	KMH4671-END-1-2	50	1.00	1.00	2.00	2.00	100.0	oilc
P4862-07	W029060-END-1-1	51	1.00	1.00	2.00	2.00	100.0	oilc
P4862-08	W029060-END-1-2	52	1.00	1.00	2.00	2.00	100.0	oilc
P4865-01	02-9042-1-1	41	1.00	1.00	2.00	2.00	100.0	oilc
P4865-02	02-9042-1-2	42	1.00	1.00	2.00	2.00	100.0	oilc
P4865-03	BC271303-1-1	43	1.00	1.00	2.00	2.00	100.0	oilc
P4865-04	BC271303-1-2	44	1.00	1.00	2.00	2.00	100.0	oilc
P4869-01	EO-02-11152024	53	1.16	8.50	9.66	8.89	90.9	
P4869-02	EO-02-11152024-E2	54	1.18	8.52	9.7	9.1	93.0	
P4870-01	TP-1	3	1.14	8.69	9.83	8.9	89.3	
P4870-02	TP-1-EPH	4	1.15	8.45	9.6	8.61	88.3	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 11/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 11/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133466

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4870-03	TP-1-VOC	5	1.12	8.75	9.87	8.84	88.2	
P4870-04	MH-735	6	1.18	8.59	9.77	8.84	89.2	
P4870-05	MH-735-EPH	7	1.18	8.42	9.6	8.68	89.1	
P4870-06	MH-735-VOC	8	1.18	8.47	9.65	8.71	88.9	
P4870-07	MH-736	9	1.19	8.63	9.82	8.85	88.8	
P4870-08	MH-736-EPH	10	1.12	8.66	9.78	8.92	90.1	
P4870-09	MH-736-VOC	11	1.16	8.80	9.96	9.1	90.2	
P4870-10	TP-15	12	1.18	8.73	9.91	8.77	86.9	
P4870-11	TP-15-EPH	13	1.19	8.71	9.9	8.73	86.6	
P4870-12	TP-15-VOC	14	1.18	8.42	9.6	8.54	87.4	
P4886-05	SVOC-GPC-BLANK	25	1.00	1.00	2.00	2.00	100.0	
P4886-06	PEST-GPC-BLANK	26	1.00	1.00	2.00	2.00	100.0	
P4886-07	PEST-GPC-BLANK-SPIKE	27	1.00	1.00	2.00	2.00	100.0	
P4886-08	PCB-GPC-BLANK	28	1.00	1.00	2.00	2.00	100.0	
P4886-09	PCB-GPC-BLANK-SPIKE	29	1.00	1.00	2.00	2.00	100.0	
P4886-10	SVOC-GPC2-BLANK	30	1.00	1.00	2.00	2.00	100.0	
P4886-11	PEST-GPC2-BLANK	31	1.00	1.00	2.00	2.00	100.0	
P4886-12	PEST-GPC2-BLANK-SPIKE	32	1.00	1.00	2.00	2.00	100.0	
P4886-13	PCB-GPC2-BLANK	33	1.00	1.00	2.00	2.00	100.0	
P4886-14	PCB-GPC2-BLANK-SPIKE	34	1.00	1.00	2.00	2.00	100.0	
P4887-01	MH-739	35	1.15	8.83	9.98	9.2	91.2	
P4887-03	MH-739-EPH	36	1.17	8.60	9.77	8.94	90.3	
P4887-04	MH-739-VOC	37	1.12	8.68	9.8	8.62	86.4	
P4887-05	MH-760	38	1.15	8.40	9.55	8.56	88.2	
P4887-07	MH-760-EPH	39	1.15	8.83	9.98	9.14	90.5	
P4887-08	MH-760-VOC	40	1.14	8.81	9.95	9.4	93.8	
P4889-01	VNJ-228	55	1.17	8.55	9.72	8.96	91.1	
P4889-02	VNJ-228-VOC	56	1.15	8.72	9.87	8.87	88.5	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/18/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 11/15/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 11/16/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133466

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
P4889-03	#72-11978	57	1.18	8.63	9.81	9.14	92.2	
P4889-04	#72-11978-VOC	58	1.13	8.76	9.89	9.41	94.5	
P4889-05	#72-11930	59	1.14	8.84	9.98	9.4	93.4	
P4889-06	#72-11930-VOC	60	1.14	8.84	9.98	9.52	94.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$