



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 11/22/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:35
In Date: 11/21/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 11/22/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133559

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P4946-01	BW-1	1	1.15	8.40	9.55	9.52	99.6	
P4946-02	BW-2	2	1.17	8.59	9.76	9.72	99.5	
P4948-01	337	3	1.15	8.82	9.97	9.53	95.0	
P4948-02	337-E2	4	1.16	8.45	9.61	9.29	96.2	
P4948-03	72-11944	5	1.15	8.44	9.59	9.35	97.2	
P4948-04	72-11944-E2	6	1.17	8.48	9.65	9.41	97.2	
P4949-03	CONTAINMENT-STONE	7	1.00	1.00	2.00	2.00	100.0	CONTAINMENT-STONE
P4949-04	CONTAINMENT-STONE-E2	8	1.00	1.00	2.00	2.00	100.0	CONTAINMENT-STONE
P4951-01	AU-05-112124	9	1.14	8.55	9.69	8.45	85.5	
P4951-02	AU-05-112124	10	1.15	8.82	9.97	9.07	89.8	
P4952-04	GAS-BUR-1614	11	1.00	1.00	2.00	2.00	100.0	oil sample
P4954-01	TR-05-112124	12	1.18	8.50	9.68	8.61	87.4	
P4954-02	TR-05-112124-E2	15	1.13	8.59	9.72	8.67	87.8	
P4954-03	TR-06-112124	13	1.16	8.43	9.59	8.4	85.9	
P4954-04	TR-06-112124-E2	14	1.14	8.82	9.96	9.1	90.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$