



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 12/19/2024

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 12/18/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:46
Out Date: 12/19/2024
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB133996

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
P5302-01	E2AR0	1	1.13	8.85	9.98	8.21	80.0	
P5302-02	E2AR1	2	1.16	8.44	9.6	8.6	88.2	
P5302-03	E2AR6	3	1.16	8.81	9.97	8.13	79.1	
P5302-04	E2AR6MS	4	1.16	8.81	9.97	8.13	79.1	
P5302-05	E2AR6MSD	5	1.16	8.81	9.97	8.13	79.1	
P5302-06	E2AR7	6	1.15	8.59	9.74	8.48	85.3	
P5302-07	E2AS0	7	1.15	8.83	9.98	8.6	84.4	
P5302-08	E2AR8	8	1.13	8.85	9.98	8.82	86.9	
P5302-09	E2AR9	9	1.17	8.55	9.72	8.79	89.1	
P5302-10	E2AS1	10	1.14	8.55	9.69	8.7	88.4	
P5302-11	VHBLK001	11	1.00	1.00	2.00	2.00	100.0	vhblk
P5302-14	E2AN7	12	1.18	8.56	9.74	8.24	82.5	
P5302-15	E2AP8	13	1.16	8.62	9.78	8.37	83.6	
P5302-16	E2AP9	14	1.14	8.42	9.56	8.57	88.2	
P5302-17	E2AQ0	15	1.15	8.40	9.55	8.83	91.4	
P5302-18	E2AQ2	16	1.13	8.72	9.85	9.33	94.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$