

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/8/2025

OVENTEMP IN Celsius(°C): 106
Time IN: 16:30
In Date: 01/07/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 01/08/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134178

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q1016-01	276	1	1.00	1.00	2.00	2.00	100.0	debris
Q1017-01	Y2309-0070-1-1	2	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-02	Y2309-0070-1-2	3	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-03	Y2309-0071-1-1	4	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-04	Y2309-0071-1-2	5	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-05	BC107266-1-1	6	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-06	BC107266-1-2	7	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-07	BC107266-2-1	8	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-08	BC107266-2-1	9	1.00	1.00	2.00	2.00	100.0	pilc
Q1017-09	HOPPER-COMP	10	1.00	1.00	2.00	2.00	100.0	pilc
Q1018-01	3732	11	1.2	8.38	9.58	8.61	88.4	
Q1018-02	3735	12	1.00	1.00	2.00	2.00	100.0	debris
Q1018-03	3744	13	1.00	1.00	2.00	2.00	100.0	debris
Q1018-04	3745	14	1.00	1.00	2.00	2.00	100.0	debris
Q1019-01	HIH072I-1-1	15	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-02	HIH072I-1-2	16	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-03	9060-1-1	17	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-04	9060-1-2	18	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-05	6-55-1-1	19	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-06	6-55-1-2	20	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-07	BC278406-1-1	21	1.00	1.00	2.00	2.00	100.0	pilc
Q1019-08	BC278406-1-2	22	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-01	BC283069-1-1	23	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-02	BC283069-1-2	24	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-03	KMH758V-1-1	25	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-04	KMH758V-1-2	26	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-05	BC270896-1-1	27	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-06	BC270896-1-2	28	1.00	1.00	2.00	2.00	100.0	pilc



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 1/8/2025

OVENTEMP IN Celsius(°C): 106
Time IN: 16:30
In Date: 01/07/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 01/08/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134178

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
Q1020-07	BC107490-1-1	29	1.00	1.00	2.00	2.00	100.0	pilc
Q1020-08	BC107490-1-2	30	1.00	1.00	2.00	2.00	100.0	pilc
Q1021-01	424	31	1.14	8.80	9.94	9.36	93.4	
Q1024-01	KMH455R-1-1	36	1.00	1.00	2.00	2.00	100.0	pilc
Q1024-02	KMH455R-1-2	37	1.00	1.00	2.00	2.00	100.0	pilc
Q1025-01	GREENBROOK-COMP	38	1.18	8.43	9.61	8.46	86.4	
Q1025-02	GREENBROOK-E2	39	1.16	8.81	9.97	8.97	88.6	
Q1027-01	EO-03-010725	32	1.15	8.74	9.89	9.2	92.1	
Q1027-02	EO-03-010725-E2	33	1.19	8.66	9.85	9.15	91.9	
Q1028-01	OR-03-010725	34	1.19	8.62	9.81	9.28	93.9	
Q1028-02	OR-03-010725-E2	35	1.12	8.70	9.82	9.22	93.1	
Q1029-01	HR-02-01072025	40	1.18	8.54	9.72	9.18	93.7	
Q1029-02	HR-02-01072025-E2	41	1.16	8.44	9.6	8.97	92.5	
Q1029-03	HR-01-01072025	42	1.19	8.65	9.84	9.21	92.7	
Q1029-04	HR-01-01072025-E2	43	1.18	8.45	9.63	9.11	93.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$