

# PERCENT SOLID

Supervisor: Iwona  
Analyst: jignesh  
Date: 1/21/2025

OVENTEMP IN Celsius(°C): 107  
Time IN: 13:50  
In Date: 01/20/2025  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:05  
Out Date: 01/21/2025  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: % SOLIDS-OVEN

QC:LB134344

| Lab ID   | Client SampleID | Dish # | Dish Wt(g) (A) | Sample Wt(g) | Dish + Sample Wt(g) (B) | Dish+Dry Sample Wt(g) (C) | % Solid | Comments |
|----------|-----------------|--------|----------------|--------------|-------------------------|---------------------------|---------|----------|
| Q1108-01 | YE8F5           | 1      | 1.13           | 8.50         | 9.63                    | 9.2                       | 94.9    |          |
| Q1108-02 | YE8F7           | 2      | 1.15           | 8.45         | 9.6                     | 8.74                      | 89.8    |          |
| Q1108-03 | YE8F7MS         | 3      | 1.15           | 8.45         | 9.6                     | 8.74                      | 89.8    |          |
| Q1108-04 | YE8F7MSD        | 4      | 1.15           | 8.45         | 9.6                     | 8.74                      | 89.8    |          |
| Q1108-05 | YE8G0           | 5      | 1.17           | 8.55         | 9.72                    | 9.13                      | 93.1    |          |
| Q1108-06 | YE8G1           | 6      | 1.15           | 8.84         | 9.99                    | 8.25                      | 80.3    |          |
| Q1108-07 | YE8G2           | 7      | 1.12           | 8.67         | 9.79                    | 9.35                      | 94.9    |          |
| Q1108-08 | YE8F6           | 8      | 1.19           | 8.72         | 9.91                    | 9.11                      | 90.8    |          |
| Q1108-09 | YE8F8           | 9      | 1.13           | 8.65         | 9.78                    | 9.18                      | 93.1    |          |
| Q1108-10 | YE8F9           | 10     | 1.17           | 8.69         | 9.86                    | 9.14                      | 91.7    |          |
| Q1108-11 | VHBLK001        | 11     | 1.00           | 1.00         | 2.00                    | 2.00                      | 100.0   | VHBLK    |
| Q1108-14 | YE8E3           | 12     | 1.18           | 8.81         | 9.99                    | 7.96                      | 77.0    |          |
| Q1108-15 | YE8E4           | 13     | 1.14           | 8.39         | 9.53                    | 7.91                      | 80.7    |          |
| Q1108-16 | YE8E5           | 14     | 1.17           | 8.67         | 9.84                    | 7.88                      | 77.4    |          |
| Q1108-17 | YE8E6           | 15     | 1.13           | 8.76         | 9.89                    | 9.4                       | 94.4    |          |
| Q1108-19 | YE8H9           | 16     | 1.14           | 8.50         | 9.64                    | 8.15                      | 82.5    |          |
| Q1108-20 | YE8J0           | 17     | 1.16           | 8.77         | 9.93                    | 8.07                      | 78.8    |          |
| Q1108-21 | YE8J1           | 18     | 1.17           | 8.69         | 9.86                    | 8.31                      | 82.2    |          |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$