

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/5/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 03/04/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/05/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134889

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q1480-01	TP-1	1	1.13	8.86	9.99	8.95	88.3	
Q1480-02	TP-1-E2	2	1.12	8.86	9.98	9.14	90.5	
Q1480-03	TP-1-E3	3	1.15	8.76	9.91	8.11	79.5	
Q1480-04	TP-1-E4	4	1.15	8.48	9.63	8.73	89.4	
Q1480-05	TP-1-E5	5	1.15	8.82	9.97	8.71	85.7	
Q1480-06	TP-1-E6	6	1.15	8.81	9.96	8.93	88.3	
Q1480-07	TP-1-E7	7	1.14	8.71	9.85	8.61	85.8	
Q1480-09	TP-2	8	1.12	8.74	9.86	8.56	85.1	
Q1480-10	TP-2-E2	9	1.14	8.82	9.96	8.53	83.8	
Q1480-11	TP-2-E3	10	1.15	8.50	9.65	8.62	87.9	
Q1480-12	TP-2-E4	11	1.14	8.76	9.9	8.8	87.4	
Q1480-13	TP-2-E5	12	1.15	8.82	9.97	8.92	88.1	
Q1480-14	TP-2-E6	13	1.16	8.82	9.98	8.95	88.3	
Q1480-15	TP-2-E7	14	1.12	8.73	9.85	8.88	88.9	
Q1480-17	TP-3	15	1.15	8.83	9.98	9.11	90.1	
Q1480-18	TP-3-E2	16	1.15	8.83	9.98	8.89	87.7	
Q1480-19	TP-3-E3	17	1.13	8.58	9.71	8.67	87.9	
Q1480-20	TP-3-E4	18	1.15	8.61	9.76	8.51	85.5	
Q1480-21	TP-3-E5	19	1.15	8.80	9.95	8.8	86.9	
Q1480-22	TP-3-E6	20	1.15	8.84	9.99	8.92	87.9	
Q1480-23	TP-3-E7	21	1.19	8.63	9.82	8.75	87.6	
Q1480-25	TP-4	22	1.14	8.83	9.97	9.07	89.8	
Q1480-26	TP-4-E2	23	1.17	8.80	9.97	9.32	92.6	
Q1480-27	TP-4-E3	24	1.13	8.76	9.89	8.95	89.3	
Q1480-28	TP-4-E4	25	1.16	8.81	9.97	9.13	90.5	
Q1480-29	TP-4-E5	26	1.12	8.67	9.79	8.97	90.5	
Q1480-30	TP-4-E6	27	1.15	8.83	9.98	9.34	92.8	
Q1480-31	TP-4-E7	28	1.11	8.62	9.73	8.88	90.1	

PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/5/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 03/04/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/05/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134889

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q1481-01	WC-K1311	29	1.15	8.70	9.85	9.00	90.2	
Q1482-01	OR-03-030425	30	1.15	8.48	9.63	8.85	90.8	
Q1482-02	OR-03-030425-E2	31	1.15	8.83	9.98	9.17	90.8	
Q1484-01	TR-OTR-COMP-01	35	1.15	8.83	9.98	9.14	90.5	
Q1484-02	TR-OTR-VOC-01	36	1.15	8.81	9.96	9.22	91.6	
Q1484-03	TR-OTR-01	37	1.12	8.63	9.75	8.83	89.3	
Q1484-04	TR-OTR-02	38	1.15	8.78	9.93	8.99	89.3	
Q1484-05	TR-OTR-03	39	1.14	8.63	9.77	8.6	86.4	
Q1485-01	DN-B-41	32	1.16	8.53	9.69	8.25	83.1	
Q1486-01	DN-B-40	33	1.16	8.40	9.56	9.02	93.6	
Q1487-01	DN-B-42	34	1.13	8.84	9.97	8.92	88.1	
Q1488-01	ENV-101-SB01	40	1.13	8.45	9.58	8.61	88.5	
Q1488-03	ENV-101-SB02	41	1.15	8.82	9.97	9.14	90.6	
Q1488-05	ENV-102-SB01	42	1.16	8.82	9.98	8.72	85.7	
Q1488-07	ENV-102-SB02	43	1.13	8.55	9.68	7.98	80.1	
Q1488-09	ENV-104-SB01	44	1.12	8.50	9.62	8.14	82.6	
Q1488-11	ENV-104-SB02	45	1.14	8.85	9.99	8.32	81.1	
Q1489-01	TP-1	46	1.13	8.86	9.99	8.18	79.6	
Q1489-02	TP-1-E2	47	1.15	8.83	9.98	8.17	79.5	
Q1489-03	TP-1-E3	48	1.13	8.70	9.83	7.68	75.3	
Q1489-05	TP-2	49	1.11	8.76	9.87	6.93	66.4	
Q1489-06	TP-2-E2	50	1.16	8.63	9.79	7.38	72.1	
Q1489-07	TP-2-E3	51	1.15	8.83	9.98	7.53	72.3	
Q1489-09	TP-3	52	1.19	8.58	9.77	7.71	76.0	
Q1489-10	TP-3-E2	53	1.15	8.84	9.99	8.7	85.4	
Q1489-11	TP-3-E3	54	1.15	8.38	9.53	7.9	80.5	
Q1489-13	TP-4	55	1.16	8.81	9.97	8.57	84.1	
Q1489-14	TP-4-E2	56	1.16	8.82	9.98	8.46	82.8	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/5/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:40
In Date: 03/04/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/05/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134889

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Sample Wt (g)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
Q1489-15	TP-4-E3	57	1.16	8.50	9.66	8.02	80.7	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$