



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/7/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 03/06/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:41
Out Date: 03/07/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134919

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q1433-03	K084-1C	1	1.12	8.76	9.88	5.68	52.1	
Q1433-06	K084-2C	2	1.15	8.48	9.63	7.29	72.4	
Q1433-09	K084-3C	3	1.15	8.80	9.95	7.89	76.6	
Q1434-02	K084-4C	4	1.15	8.62	9.77	7.38	72.3	
Q1434-05	K084-5C	5	1.16	8.81	9.97	8.68	85.4	
Q1434-07	K084-6C	6	1.12	8.70	9.82	7.49	73.2	
Q1434-10	K084-7C	7	1.12	8.76	9.88	7.23	69.7	
Q1436-02	K084-8B	8	1.17	8.80	9.97	8.01	77.7	
Q1436-03	K084-8C	9	1.16	8.82	9.98	7.63	73.4	
Q1436-05	K084-9B	10	1.16	8.82	9.98	7.12	67.6	
Q1436-06	K084-9C	11	1.16	8.67	9.83	6.26	58.8	
Q1436-08	K084-10B	12	1.16	8.57	9.73	7.16	70.0	
Q1436-09	K084-10C	13	1.16	8.80	9.96	6.66	62.5	
Q1437-01	K084-11B	14	1.13	8.80	9.93	7.54	72.8	
Q1437-02	K084-11C	15	1.15	8.82	9.97	6.59	61.7	
Q1437-04	K084-12B	16	1.13	8.74	9.87	7.61	74.1	
Q1437-05	K084-12C	17	1.15	8.82	9.97	7.03	66.7	
Q1437-07	K084-13B	18	1.19	8.72	9.91	7.45	71.8	
Q1437-08	K084-13C	19	1.16	8.45	9.61	6.95	68.5	
Q1437-10	K084-14B	20	1.16	8.80	9.96	7.03	66.7	
Q1438-01	K084-14C	21	1.15	8.62	9.77	7.03	68.2	
Q1438-03	K084-15C	22	1.14	8.80	9.94	7.71	74.7	
Q1438-05	K084-16C	23	1.19	8.53	9.72	7.56	74.7	
Q1438-09	K084-DUP4	24	1.18	8.81	9.99	6.93	65.3	
Q1507-01	50-MIDDLESEX-AVE	25	1.15	8.83	9.98	8.92	88.0	
Q1507-02	50-MIDDLESEX-AVE-E2	26	1.15	8.80	9.95	8.41	82.5	
Q1508-01	RBR251372	27	1.12	8.87	9.99	7.42	71.0	
Q1508-02	RBR251372-E2	28	1.19	8.75	9.94	7.76	75.1	



PERCENT SOLID

Supervisor: Iwona
Analyst: jignesh
Date: 3/7/2025

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 03/06/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN#1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:41
Out Date: 03/07/2025
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLID- OVEN

QC:LB134919

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Sample Wt(g)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
Q1509-01	HIH088X-1-1	29	1.00	1.00	2.00	2.00	100.0	oilc
Q1509-02	HIH088X-1-2	30	1.00	1.00	2.00	2.00	100.0	oilc
Q1510-01	FMC-25-0001-0005	31	1.16	8.50	9.66	7.99	80.4	
Q1510-02	FMC-25-0006	32	1.00	1.00	2.00	2.00	100.0	oil sample
Q1514-01	ENV-105-SB01	39	1.15	8.84	9.99	9.25	91.6	
Q1514-03	ENV-105-SB02	40	1.13	8.84	9.97	8.29	81.0	
Q1514-05	ENV-103-SB01	41	1.14	8.81	9.95	8.99	89.1	
Q1515-01	AU-06-030625	42	1.15	8.81	9.96	8.88	87.7	
Q1515-02	AU-06-030625-E2	43	1.16	8.82	9.98	8.62	84.6	
Q1516-01	GAS-AUD-25-0004	44	1.17	8.58	9.75	9.24	94.1	
Q1516-02	GAS-AUD-25-0025	45	1.16	8.68	9.84	9.7	98.4	
Q1518-01	OILY-DEBRIS-COMP	46	1.15	8.83	9.98	9.29	92.2	
Q1520-01	BC248121-1-1	33	1.00	1.00	2.00	2.00	100.0	oilc
Q1520-02	BC248121-1-2	34	1.00	1.00	2.00	2.00	100.0	oilc
Q1520-03	BC260492-1-1	35	1.00	1.00	2.00	2.00	100.0	oilc
Q1520-04	BC260492-1-2	36	1.00	1.00	2.00	2.00	100.0	oilc
Q1520-05	JEH488H-1-1	37	1.00	1.00	2.00	2.00	100.0	oilc
Q1520-06	JEH488H-1-2	38	1.00	1.00	2.00	2.00	100.0	oilc

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$