


Ident	Instrument IC-1				Analyst: NF		Method: 300.0 / 9056A		Initial Analyst	
	Con F-	Con CL-	Con NO2	Con BR-	Con NO3	Con HPO4	Con SO4	Method name		date time
STD1	0.000	0.000	0.000	0.000	0.000	0.000	0.000	IC1-042225	4/22/2025 10:37	10 NF/IZ
STD2	0.431	0.651	0.655	2.172	0.551	1.090	3.366	IC1-042225	4/22/2025 10:59	10 NF/IZ
STD3	0.810	1.217	1.222	4.049	1.011	2.010	6.091	IC1-042225	4/22/2025 11:20	10 NF/IZ
STD4	0.984	1.509	1.498	5.011	1.250	2.445	7.369	IC1-042225	4/22/2025 11:42	10 NF/IZ
STD5	1.938	2.914	2.916	9.735	2.424	4.860	14.361	IC1-042225	4/22/2025 12:03	10 NF/IZ
STD6	4.065	5.918	5.915	19.752	4.949	10.183	30.380	IC1-042225	4/22/2025 12:25	10 NF/IZ
STD7	4.972	7.591	7.594	25.280	6.316	12.412	36.933	IC1-042225	4/22/2025 12:46	10 NF/IZ
ICV	2.034	3.100	3.074	10.246	2.611	5.165	15.120	IC1-042225	4/22/2025 13:07	10 NF/IZ
ICB	0.000	0.000	0.000	0.000	0.000	0.000	0.000	IC1-042225	4/22/2025 13:58	10 NF/IZ
CCV	2.062	3.094	3.102	10.398	2.568	5.056	15.119	IC1-042225	5/21/2025 10:38	10 NF/IZ
CCB	0.000	0.000	0.000	0.000	0.000	0.000	0.000	IC1-042225	5/21/2025 11:00	10 NF/IZ
LB135856BLW	0.000	0.000	0.000	0.000	0.000	0.000	0.000	IC1-042225	5/21/2025 11:21	10 NF/IZ
LB135856BSW	2.053	3.087	3.113	10.382	2.564	5.099	15.069	IC1-042225	5/21/2025 11:43	10 NF/IZ
Q2088-01	0.189	9.967	0.096	0.000	2.212	0.000	42.443	IC1-042225	5/21/2025 12:04	10 NF/IZ
Q2093-14	0.188	94.984	0.179	0.426	0.448	0.000	61.086	IC1-042225	5/21/2025 12:26	10 NF/IZ
Q2093-14MS	2.138	93.977	3.099	10.074	2.874	5.020	74.590	IC1-042225	5/21/2025 12:47	10 NF/IZ
Q2093-14MSD	2.209	94.156	3.237	10.531	2.985	5.051	75.303	IC1-042225	5/21/2025 13:09	10 NF/IZ
Q2102-03	0.124	9.211	0.000	0.000	0.000	0.000	6.882	IC1-042225	5/21/2025 13:30	10 NF/IZ
CCV	2.076	3.097	3.113	10.386	2.570	4.963	15.130	IC1-042225	5/21/2025 14:07	10 NF/IZ
CCB	0.000	0.000	0.000	0.000	0.000	0.000	0.000	IC1-042225	5/21/2025 14:29	10 NF/IZ

Clear table

Instrument ID: IC-2 Analyst: IZ Method: 300.0 / 9056A									
ident	concentra tion F-	concentratio n CL-	concentrati on NO2	concentrati on BR-	concentrati on NO3	concentrati on HPO4	concentrati on SO4	file name	Initial wt/ Final
STD1	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	IC1-042225	10 NF/IZ
STD2	0.4310	0.6510	0.6550	2.1720	0.5510	1.0900	3.3660	IC1-042225	10 NF/IZ
STD3	0.8100	1.2170	1.2220	4.0490	1.0110	2.0100	6.0910	IC1-042225	10 NF/IZ
STD4	0.9840	1.5090	1.4980	5.0110	1.2500	2.4450	7.3690	IC1-042225	10 NF/IZ
STD5	1.9380	2.9140	2.9160	9.7350	2.4240	4.8600	14.3610	IC1-042225	10 NF/IZ
STD6	4.0650	5.9180	5.9150	19.7520	4.9490	10.1830	30.3800	IC1-042225	10 NF/IZ
STD7	4.9720	7.5910	7.5940	25.2800	6.3160	12.4120	36.9330	IC1-042225	10 NF/IZ
ident	True Value F-	True Value CL-	True Value NO2	True Value BR-	True Value NO3	True Value HPO4	True Value SO4		
STD1	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000		
STD2	0.4000	0.6000	0.6000	2.0000	0.5000	1.0000	3.0000		
STD3	0.8000	1.2000	1.2000	4.0000	1.0000	2.0000	6.0000		
STD4	1.0000	1.5000	1.5000	5.0000	1.2500	2.5000	7.5000		
STD5	2.0000	3.0000	3.0000	10.0000	2.5000	5.0000	15.0000		
STD6	4.0000	6.0000	6.0000	20.0000	5.0000	10.0000	30.0000		
STD7	5.0000	7.5000	7.5000	25.0000	6.2500	12.5000	37.0000		
ident	Relative Error F-	Relative Error CL-	Relative Error NO2	Relative Error BR-	Relative Error NO3	Relative Error HPO4	Relative Error SO4		
STD1									
STD2	7.7500	8.5000	9.1667	8.6000	10.2000	9.0000	12.2000		
STD3	1.2500	1.4167	1.8333	1.2250	1.1000	0.5000	1.5167		
STD4	-1.6000	0.6000	-0.1333	0.2200	0.0000	-2.2000	-1.7467		
STD5	-3.1000	-2.8667	-2.8000	-2.6500	-3.0400	-2.8000	-4.2600		
STD6	1.6250	-1.3667	-1.4167	-1.2400	-1.0200	1.8300	1.2667		
STD7	-0.5600	1.2133	1.2533	1.1200	1.0560	-0.7040	-0.1811		

Fluoride (Anions)




Function: $A = -6.13443\text{E-}3 + 0.0144025 \times Q$

Relative standard deviation 2.337762 %

Correlation coefficient 0.999714

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	0.400	10.0	1.0	1.0	0.056	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	0.800	10.0	1.0	1.0	0.110	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	1.000	10.0	1.0	1.0	0.136	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	2.000	10.0	1.0	1.0	0.273	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	4.000	10.0	1.0	1.0	0.579	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	5.000	10.0	1.0	1.0	0.710	STD7	2025-04-22 12:46:28 UTC-4	used

Chloride (Anions)



Function: $A = -5.12446\text{E-}3 + 9.72422\text{E-}3 \times Q$


Relative standard deviation 2.452549 %

Correlation coefficient 0.999683

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	0.600	10.0	1.0	1.0	0.058	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	1.200	10.0	1.0	1.0	0.113	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	1.500	10.0	1.0	1.0	0.142	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	3.000	10.0	1.0	1.0	0.278	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	6.000	10.0	1.0	1.0	0.570	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	7.500	10.0	1.0	1.0	0.733	STD7	2025-04-22 12:46:28 UTC-4	used

Nitrite (Anions)

($\mu\text{S}/\text{cm}$) $\times$ min



Function: $A = -0.0187238 + 0.0219021 \times Q$


Relative standard deviation 2.531721 %

Correlation coefficient 0.999669

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	0.600	10.0	1.0	1.0	0.125	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	1.200	10.0	1.0	1.0	0.249	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	1.500	10.0	1.0	1.0	0.309	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	3.000	10.0	1.0	1.0	0.620	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	6.000	10.0	1.0	1.0	1.277	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	7.500	10.0	1.0	1.0	1.644	STD7	2025-04-22 12:46:28 UTC-4	used

Bromide (Anions)

($\mu\text{S}/\text{cm}$) $\times$ min



Function: $A = -9.36439E-3 + 4.26318E-3 \times Q$


Relative standard deviation: 2.280508 %

Correlation coefficient: 0.999728

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	2.000	10.0	1.0	1.0	0.083	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	4.000	10.0	1.0	1.0	0.163	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	5.000	10.0	1.0	1.0	0.204	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	10.000	10.0	1.0	1.0	0.406	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	20.000	10.0	1.0	1.0	0.833	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	25.000	10.0	1.0	1.0	1.068	STD7	2025-04-22 12:46:28 UTC-4	used

Nitrate (Anions)

($\mu\text{S}/\text{cm}$) $\times$ min



Function: $A = -0.0163950 + 0.0241967 \times Q$


Relative standard deviation: 2.309177 %

Correlation coefficient 0.999724

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	0.500	10.0	1.0	1.0	0.117	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	1.000	10.0	1.0	1.0	0.228	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	1.250	10.0	1.0	1.0	0.286	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	2.500	10.0	1.0	1.0	0.570	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	5.000	10.0	1.0	1.0	1.181	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	6.250	10.0	1.0	1.0	1.512	STD7	2025-04-22 12:46:28 UTC-4	used

Phosphate (Anions)


(μS/cm) x min

Function: $A = -9.96993E-3 + 8.43655E-3 \times Q$

Relative standard deviation 2.492252 %

Correlation coefficient 0.999676

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	1.000	10.0	1.0	1.0	0.082	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	2.000	10.0	1.0	1.0	0.160	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	2.500	10.0	1.0	1.0	0.196	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	5.000	10.0	1.0	1.0	0.400	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	10.000	10.0	1.0	1.0	0.849	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	12.500	10.0	1.0	1.0	1.037	STD7	2025-04-22 12:46:28 UTC-4	used

Sulfate (Anions) $(\mu\text{S}/\text{cm}) \times \text{min}$ Function: $A = -0.0453007 + 7.61540E-3 \times Q$

Relative standard deviation 2.675209 %

Correlation coefficient 0.999634


Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2025-04-22 10:37:58 UTC-4	used
Standard 2	1	3.000	10.0	1.0	1.0	0.211	STD2	2025-04-22 10:59:20 UTC-4	used
Standard 3	1	6.000	10.0	1.0	1.0	0.419	STD3	2025-04-22 11:20:44 UTC-4	used
Standard 4	1	7.500	10.0	1.0	1.0	0.516	STD4	2025-04-22 11:42:09 UTC-4	used
Standard 5	1	15.000	10.0	1.0	1.0	1.048	STD5	2025-04-22 12:03:34 UTC-4	used
Standard 6	1	30.000	10.0	1.0	1.0	2.268	STD6	2025-04-22 12:25:01 UTC-4	used
Standard 7	1	37.000	10.0	1.0	1.0	2.767	STD7	2025-04-22 12:46:28 UTC-4	used

Sample data

Ident STD1
Sample type Standard 1
Determination start 2025-04-22 10:37:58 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.32 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions


Sample data

Ident STD2
Sample type Standard 2
Determination start 2025-04-22 10:59:20 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.26 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.363	0.0560	0.425	0.431	Fluoride
2	6.502	0.0581	0.371	0.651	Chloride
3	8.570	0.1247	0.628	0.655	Nitrite
4	10.562	0.0832	0.367	2.172	Bromide
5	12.222	0.1168	0.452	0.551	Nitrate
6	14.197	0.0820	0.201	1.090	Phosphate
7	15.773	0.2111	0.575	3.366	Sulfate


Sample data

Ident STD3
Sample type Standard 3
Determination start 2025-04-22 11:20:44 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.43 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.360	0.1105	0.845	0.810	Fluoride
2	6.503	0.1132	0.722	1.217	Chloride
3	8.577	0.2488	1.250	1.222	Nitrite
4	10.572	0.1632	0.724	4.049	Bromide
5	12.233	0.2282	0.886	1.011	Nitrate
6	14.192	0.1596	0.393	2.010	Phosphate
7	15.757	0.4185	1.139	6.091	Sulfate


Sample data

Ident STD4
Sample type Standard 4
Determination start 2025-04-22 11:42:09 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.43 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.355	0.1356	1.044	0.984	Fluoride
2	6.502	0.1416	0.903	1.509	Chloride
3	8.577	0.3094	1.559	1.498	Nitrite
4	10.573	0.2043	0.902	5.011	Bromide
5	12.235	0.2860	1.108	1.250	Nitrate
6	14.177	0.1963	0.510	2.445	Phosphate
7	15.743	0.5159	1.419	7.369	Sulfate


Sample data

Ident STD5
 Sample type Standard 5
 Determination start 2025-04-22 12:03:34 UTC-4
 Method IC1-042225
 Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
 Channel Conductivity
 Recording time 19.0 min
 Integration Automatically
 Column type Metrosep A Supp 19 - 150/4.0
 Eluent composition not defined
 Flow 0.700 mL/min
 Maximum flow monitored yes
 Pressure 11.49 MPa
 Maximum pressure monitored yes
 Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.360	0.2730	2.125	1.938	Fluoride
2	6.502	0.2783	1.784	2.914	Chloride
3	8.573	0.6200	3.144	2.916	Nitrite
4	10.562	0.4057	1.810	9.735	Bromide
5	12.218	0.5702	2.226	2.424	Nitrate
6	14.170	0.4000	1.020	4.860	Phosphate
7	15.748	1.0484	2.866	14.361	Sulfate


Sample data

Ident STD6
Sample type Standard 6
Determination start 2025-04-22 12:25:01 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.49 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.360	0.5793	4.564	4.065	Fluoride
2	6.497	0.5704	3.686	5.918	Chloride
3	8.570	1.2768	6.510	5.915	Nitrite
4	10.550	0.8327	3.746	19.752	Bromide
5	12.202	1.1810	4.667	4.949	Nitrate
6	14.167	0.8491	2.199	10.183	Phosphate
7	15.753	2.2683	6.294	30.380	Sulfate


Sample data

Ident STD7
Sample type Standard 7
Determination start 2025-04-22 12:46:28 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.43 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.353	0.7100	5.662	4.972	Fluoride
2	6.480	0.7331	4.785	7.591	Chloride
3	8.547	1.6445	8.440	7.594	Nitrite
4	10.517	1.0684	4.850	25.280	Bromide
5	12.162	1.5118	6.035	6.316	Nitrate
6	14.127	1.0372	2.803	12.412	Phosphate
7	15.722	2.7673	7.822	36.933	Sulfate


Sample data

Ident ICV
Sample type Check standard 1
Determination start 2025-04-22 13:07:56 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.60 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.362	0.2867	2.234	2.034	Fluoride
2	6.493	0.2963	1.906	3.100	Chloride
3	8.555	0.6546	3.328	3.074	Nitrite
4	10.532	0.4274	1.914	10.246	Bromide
5	12.185	0.6154	2.414	2.611	Nitrate
6	14.168	0.4257	1.100	5.165	Phosphate
7	15.775	1.1062	3.034	15.120	Sulfate


Sample data

Ident ICB
Sample type Sample
Determination start 2025-04-22 13:58:31 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.32 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time	Area	Height	Concentration	Component name
	min	(µS/cm) x min	µS/cm	ppm	
1	1.362	0.0024	0.015	invalid	


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2025-05-21 10:38:35 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.50 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.348	0.2908	2.242	2.062	Fluoride
2	6.585	0.2957	1.860	3.094	Chloride
3	8.735	0.6607	3.247	3.102	Nitrite
4	10.832	0.4339	1.864	10.398	Bromide
5	12.512	0.6049	2.286	2.568	Nitrate
6	13.685	0.4166	1.096	5.056	Phosphate
7	15.493	1.1061	3.032	15.119	Sulfate


Sample data

Ident CCB
Sample type Sample
Determination start 2025-05-21 11:00:05 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.50 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Sample data

Ident LB135856BLW
Sample type Sample
Determination start 2025-05-21 11:21:35 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.39 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions


Sample data

Ident LB135856BSW
Sample type Check standard 1
Determination start 2025-05-21 11:43:05 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.84 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.347	0.2896	2.246	2.053	Fluoride
2	6.565	0.2951	1.866	3.087	Chloride
3	8.705	0.6631	3.272	3.113	Nitrite
4	10.787	0.4332	1.872	10.382	Bromide
5	12.463	0.6039	2.293	2.564	Nitrate
6	13.670	0.4202	1.113	5.099	Phosphate
7	15.492	1.1023	3.028	15.069	Sulfate


Sample data

Ident Q2088-01
Sample type Sample
Determination start 2025-05-21 12:04:37 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.56 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.345	0.0210	0.149	0.189	Fluoride
2	6.573	0.9641	6.185	9.967	Chloride
3	7.345	0.0109	0.024	invalid	
4	8.702	0.0024	0.012	0.096	Nitrite
5	12.468	0.5189	1.972	2.212	Nitrate
6	15.493	3.1869	8.936	42.443	Sulfate


Sample data

Ident Q2093-14
Sample type Sample
Determination start 2025-05-21 12:26:09 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.39 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.350	0.0209	0.148	0.188	Fluoride
2	6.595	9.2313	62.212	94.984	Chloride
3	7.308	0.0060	0.017	invalid	
4	8.708	0.0206	0.101	0.179	Nitrite
5	10.202	0.0088	0.037	0.426	Bromide
6	12.470	0.0921	0.347	0.448	Nitrate
7	15.492	4.6066	13.028	61.086	Sulfate


Sample data

Ident Q2093-14MS
Sample type Sample
Determination start 2025-05-21 12:47:42 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.61 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.342	0.3018	2.343	2.138	Fluoride
2	6.580	9.1334	61.602	93.977	Chloride
3	7.277	0.0061	0.016	invalid	
4	8.693	0.6600	3.292	3.099	Nitrite
5	10.173	0.0074	0.034	invalid	
6	10.768	0.4201	1.835	10.074	Bromide
7	12.442	0.6790	2.593	2.874	Nitrate
8	13.650	0.4135	1.092	5.020	Phosphate
9	15.455	5.6351	16.046	74.590	Sulfate


Sample data

Ident Q2093-14MSD
Sample type Sample
Determination start 2025-05-21 13:09:16 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.39 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.343	0.3120	2.424	2.209	Fluoride
2	6.583	9.1508	61.756	94.156	Chloride
3	7.278	0.0060	0.016	invalid	
4	8.700	0.6904	3.439	3.237	Nitrite
5	10.190	0.0071	0.033	invalid	
6	10.777	0.4396	1.921	10.531	Bromide
7	12.450	0.7058	2.703	2.985	Nitrate
8	13.660	0.4162	1.092	5.051	Phosphate
9	15.468	5.6893	16.212	75.303	Sulfate


Sample data

Ident Q2102-03
Sample type Sample
Determination start 2025-05-21 13:30:51 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.33 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	2.910	0.0069	0.030	invalid	
2	3.352	0.0117	0.081	0.124	Fluoride
3	3.573	0.0017	0.009	invalid	
4	6.573	0.8906	5.721	9.211	Chloride
5	7.377	0.0047	0.011	invalid	
6	15.495	0.4788	1.303	6.882	Sulfate


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2025-05-21 14:07:57 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.50 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.352	0.2929	2.236	2.076	Fluoride
2	6.590	0.2961	1.872	3.097	Chloride
3	8.742	0.6630	3.268	3.113	Nitrite
4	10.840	0.4334	1.871	10.386	Bromide
5	12.523	0.6054	2.294	2.570	Nitrate
6	13.690	0.4088	1.074	4.963	Phosphate
7	15.498	1.1069	3.039	15.130	Sulfate

Sample data

Ident CCB
Sample type Sample
Determination start 2025-05-21 14:29:16 UTC-4
Method IC1-042225
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.56 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions