

Analysis Method: 9034

ANALYST: rubina

Parameter: Reactive Sulfide

SUPERVISOR REVIEW BY: Iwona

Run Number: LB135891

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3114
Starch Solution, 4L	W3149

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
1	PB168111BL		1	5.00	50	2.00	0.00	1.92	1.92	0.08	0.00	0.00	05/22/2025	16:15
2	Q2092-02		1	5.02	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	05/22/2025	16:17
3	Q2092-02DUP		1	5.02	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	05/22/2025	16:19
4	Q2093-02		1	5.03	50	2.00	0.00	1.86	1.86	0.14	0.06	4.77	05/22/2025	16:21
5	Q2095-04		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	05/22/2025	16:23
6	Q2095-08		1	5.04	50	2.00	0.00	1.88	1.88	0.12	0.04	3.17	05/22/2025	16:25
7	Q2097-02		1	5.02	50	2.00	0.00	1.86	1.86	0.14	0.06	4.78	05/22/2025	16:28
8	Q2097-04		1	5.07	50	2.00	0.00	1.88	1.88	0.12	0.04	3.16	05/22/2025	16:30
9	Q2097-06		1	5.05	50	2.00	0.00	1.88	1.88	0.12	0.04	3.17	05/22/2025	16:33
10	Q2097-08		1	5.01	50	2.00	0.00	1.86	1.86	0.14	0.06	4.79	05/22/2025	16:36
11	Q2097-10		1	5.01	50	2.00	0.00	1.86	1.86	0.14	0.06	4.79	05/22/2025	16:38
12	Q2097-12		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	05/22/2025	16:40

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume

Analysis Method: 9034

Parameter: Reactive Sulfide

Run Number: LB135891

ANALYST: rubina

SUPERVISOR REVIEW BY: Iwona

Constant: 16000

Normality1: 0.025

Normality2: 0.025

Reagent/Standard	Lot/Log #
SODIUM THIOSULFATE, 0.025N, 4LITRE	W3105
IODINE SOLUTION .025N 1L	W3114
Starch Solution, 4L	W3149

Seq	Lab ID	True Value (mg/l)	DF	Initial Weight (g)	Final Volume (ml)	T1 (ml)	T2 Initial	T2 Final	T2 Diff. (ml)	T1 - T2 Diff (mL)	Value Corrected With Blank	Result (ppm)	Anal Date	Anal Time
13	Q2097-14		1	5.04	50	2.00	0.00	1.88	1.88	0.12	0.04	3.17	05/22/2025	16:43
14	Q2097-16		1	5.08	50	2.00	0.00	1.86	1.86	0.14	0.06	4.72	05/22/2025	16:46
15	Q2097-18		1	5.04	50	2.00	0.00	1.90	1.90	0.10	0.02	1.59	05/22/2025	16:48
16	Q2099-01		1	5.05	50	2.00	0.00	1.90	1.90	0.10	0.02	1.58	05/22/2025	16:50
17	Q2101-04		1	5.06	50	2.00	0.00	1.86	1.86	0.14	0.06	4.74	05/22/2025	16:52
18	Q2109-04		1	5.01	50	2.00	0.00	1.90	1.90	0.10	0.02	1.60	05/22/2025	16:53

T1 = Titrant1

T2 = Titrant2

T2 Diff = T2 Final - T2 Initial

Value Corrected With Blank = ((T1 - T2 Diff) - Blank Correction(BL))

Result = ((T1 * Normality1) - ((T1 - Value Corrected With Blank) * Normality2)) * Constant / Initial Volume